**Mode d'emploi
PIECES A MAIN
PM2**

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**Operating
Instructions PM2
HANDPIECES**

Please read carefully and
keep this copy for further use

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
PM2 HANDPIECES (EN)

1. TYPE

Straight and angled handpieces with external or internal spray, without light for microsurgery.

2. SYMBOLS USED

	CE Marking with the number of the notified body.
Rx Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	CAUTION! Refer to the accompanying documents
	Materials to be recycled The disposal and/or recycling of materials must be performed in accordance with the directives and the legislation in force.
	Sterilizable in autoclave up to the specified temperature
	Cleaning in machine authorized
	Manufacturer
	Lock handpiece
	Handpiece opening direction

3. INTENDED USE

PM2 handpieces are intended for cutting and shaping bones including the spine and cranium.

4. INDICATION FOR USE

Product intended for professional use. PM2 handpieces are indicated for the drilling of bone in the head & neck/ENT (Otolgic, Neurologic, Neurotologic, Sinus, Rhinologic) Neurosurgical (cranial), and spinal surgical procedures.

5. CONTRAINDICATION

None currently known.

6. PRECAUTIONS - WARNINGS

6.1. Warnings, Precautions for use

Technical specifications, illustrations and dimensions contained in this document are for information purposes only. This information shall not be used in any claims. For small improvements, the manufacturer will not update the present instructions. For major ones, these instructions must be updated. If the regulations or standards vary from one country or region to another, the stricter rule shall prevail and should be implemented by the manufacturer.

For additional information, please contact Bien-Air Surgery SA at the address indicated on the back-cover of this document.

The device and its accessories should be used only by duly trained and competent medical personnel, in particular in compliance with the legal provisions in force regarding occupational safety, health and accident prevention measures, and the present user manual. According to these measures, the user has the following obligations:

- To only use devices that are in perfect working condition. In the event of irregular operation, excessive vibrations, abnormal overheating or other signs suggesting malfunctioning of the device, work must be suspended immediately. In this case, contact a repair center approved by Bien-Air Surgery.
- Make sure that the device is used only for the purpose for which it is intended. Protect

yourself, patients and third parties from all danger and avoid contamination by the product.

The device and its accessories are designed solely for medical treatment. Any use that does not comply with its intended use is unauthorized and may prove dangerous. This medical device complies with European legal provisions in force.

Install the handpiece on an appropriate mount to prevent risk of injury or infection to you, the patient or third parties. Use only Bien-Air Surgery SA original maintenance products, accessories and/or spare parts. The use of other products, accessories or parts could void the guarantee and/or endanger the patient or the operator.

CAUTION

Use proper irrigation fluid. Irrigants with preservatives are contraindicated for use on nervous tissue, specialists must verify their adequation before use. This is the hospital's responsibility.

6.2. Environmental protection and indications for device disposal

 This equipment must be recycled. The disposal and/or recycling of materials must be performed in accordance with the directives and the legislation in force. The user can return the device to their distributor or directly call a firm certified to treat and recover this type of equipment (European directive 2012/19/EU).

7. TECHNICAL DATA

Subject to modification of the models and technical amendments.

Environmental conditions:

Work	Temperature:	+10°C to +30°C (+50°F to +86°F)
	Relative humidity:	20% to 80%, including condensation
	Atmospheric pressure:	700 hPa to 1060 hPa

Transport	Temperature:	-25°C to +70°C (-13°F to +158°F)
	Relative humidity: Atmospheric pressure:	10% to 100%, including condensation 600 hPa to 1060 hPa
Storage	Temperature:	0°C to +40°C (+32°F to +104°F)
	Relative humidity: Atmospheric pressure:	10% to 90%, including condensation 600 hPa to 1060 hPa

Type: Straight and angled handpieces without light for microsurgery.

Rotation speed: PM2 handpieces have a ratio of 1:1.

Noise level: As per ISO 11498: < 72 dBA at 45 cm.

Size / Weight:

Handpiece reference	Designation	Length – handpiece without rotary instrument (mm)	Hand-piece weight (g)
REF 1600831	PM2-S70HD	74	40
REF 1600830	PM2-A70HD	86	49
REF 1600767	PM2-S70	74	38
REF 1600977	PM2-S125	129	45
REF 1600978	PM2-S150	154	48
REF 1600765	PM2-A70	86	46
REF 1600768	PM2-A95	111	49
REF 1600769	PM2-A125	141	53
REF 1600766	PM2-A150	166	56
REF 1600880	Cranio-guard 10mm rot.	61	23
REF 1600879	Cranio-guard 15mm rot.	66	23
REF 1600912	Cranio-guard 10mm fix	61	23
REF 1600913	Cranio-guard 15mm fix	66	23
REF 1600887	Cranio-guard 25mm fix	76	24

Period of operation: See Chap 8.8.

Classification: IIa as per 93/42 EEC directive.

Sterilization: See Chap 9.6.

Cranio-guard identification: The marking according to **Fig. 1** helps to clearly see the cranio-guard length and corresponding rotary instrument:

- Length 10mm: No line on the rotary instrument and cranio-guard.
- Length 15mm: One line on the rotary instrument and cranio-guard (**Fig. 1**).
- Length 25mm: Two lines on the rotary instrument and the cranio-guard.

For further information concerning microsurgery instruments and their accessories, please contact your Bien-Air Surgery SA local distributor or consult our homepage www.bienair.com.

8. USE

8.1. Type of coupling

Bien-Air Surgery coupling.

8.2. Drive

By **NANO**, **RAPIDO** or **BASCH-1** Bien-Air Surgery micromotor.

8.3. Type of rotary instrument

Bien-Air Surgery rotary instruments from the PM2 range must only be used with Bien-Air Surgery handpieces from the PM2 range (**Fig. 2**). Bien-Air Surgery fully declines all liability if rotary instruments made by another brand are used.

8.4. Insertion and removal of a rotary instrument

 **CAUTION**

Never insert a tool on a rotating instrument.

For all PM2 except cranio-guard and 70 “Heavy Duty”:

1. To open the handpiece, pull the ring up to the thrust stop and turn it in the direction of the arrow (**Fig. 3**).

2. Insert the rotary instrument up to the thrust stop (**Fig. 4A**).
3. To change the rotary instrument exposure to the second or third position, pull the rotary instrument and align the engraved ring with the handpiece nose.
4. To close the handpiece, turn the ring in the direction of the arrow, up to the thrust stop and release it. The ring automatically goes back to the front of the handpiece (**Fig 4B**).

PM2 80K rotary instruments can be used in every exposition up to 80,000 rpm, but for optimum drill precision, we advise that the rotary instrument is fully inserted in the handpiece at the minimum exposure, when used at speeds above 50,000 rpm.

For “Heavy Duty” PM2-S70HD and PM2-A70HD (REF 1600831 and 1600830)

1. To open the handpiece, pull the ring up to the thrust stop and turn it in the direction of the arrow (**Fig. 3**).
2. If necessary, remove the cranio-guard nose by pulling it (**Fig. 5A**).
3. Insert the nose up to the thrust stop (**Fig. 5B**).
4. Insert the rotary instrument up to the thrust stop (**Fig. 4A**).
5. To change the rotary instrument exposure to the second or third position, pull the rotary instrument and align the engraved ring with the handpiece nose.
6. To close the handpiece, turn the ring in the direction of the arrow up to the thrust stop and release it. The ring automatically returns to the front of the handpiece (**Fig. 4B**).

For cranio-guard:

 **CAUTION**

Cranio-guard must only be used with straight 70mm « Heavy Duty » handpieces (REF 1600831).

1. To open the handpiece, pull the ring up to the thrust stop and turn it in the direction of the arrow (**Fig. 3**).
2. Remove the « Heavy Duty » handpiece or cranio-guard nose by pulling it (**Fig. 5**).
3. Insert the cranio-guard rotary instrument up to the thrust stop (**Fig. 6A**).
4. Insert the cranio-guard nose onto the « Heavy Duty » handpiece (**Fig. 6B**). To fix the cranio-guard, align the mark with the irrigation channel of the « Heavy Duty » handpiece (**Fig 6C**).
5. To close the handpiece, turn the ring in the direction of the arrow, up to the thrust stop and release it. The ring automatically goes back to the front of the handpiece (**Fig. 6D**).

If the handpiece closing ring is open or not closed properly, the handpiece will automatically stop rotating when the motor is turned on.

8.5. Check

Prior to start-up, make sure that the instrument is secure by pulling it (**Fig. 4C**). Prior to working on a patient, pull the PM2 « Heavy Duty » and cranio-guard nose to be sure that it is secure, then run the rotary instrument to check that there is no vibration.

CAUTION

The rotary instrument must be in perfect condition: Shape, surface condition and cut. The user is fully responsible for:

- Ensuring that the rotary instrument meets the Bien-Air Surgery SA specifications.
- Ensuring that it rotates without vibration at the speed at which it will be used.

8.6. « Heavy Duty » handpiece irrigation

PM2-S70HD and PM2-A70HD handpiece has internal irrigation. Connect the irrigation line on the top of the handpiece (**Fig. 7A**). Irrigation can be orientated at 360° by turning the handpiece nose (**Fig. 7B**).

8.7. Precautions – warnings

Do not exceed the maximum rotation speed according to the rotary instrument range and its exposition, as stated in the PM2 rotary instrument instructions for use.

Never mount an instrument on a rotating motor.

Never adjust the clamping/unclamping ring while the instrument is rotating, since this could damage the mechanism or destroy the motor.

CAUTION

Before using the instrument, turn it on at a moderate speed with a tool in the rotary instrument locking mechanism for a few seconds, so as to spread the lubricant and remove any excess.

8.8. Operating period

To avoid exceeding 41°C on the outer surface of the handpiece when used at full power and at maximum speed, the following rules should be complied with:

- a) Limit the maximum speed of rotation based on the intended use.
- b) Adequate irrigation is recommended.
- c) Period of continuous operation highly depends on irrigation and the load and speed applied to the tool.
Indicative values:
10-15 minutes at low load.
1-2 minutes at high load without irrigation followed by a minimum interruption period of 30 seconds.

CAUTION

If the irrigation tube becomes clogged during a surgical operation, use a cleaning wire or a long cleaning curette to clear it.

9. CLEANING / MAINTENANCE / STERILIZATION

CAUTION

These devices are delivered "non-sterile".

Before use please comply with these instructions.

Precautions for use

Hospital procedures must be followed. Hospital personnel working with contaminated or potentially contaminated medical instruments must comply with universal precautions. Pointed and sharp instruments should be handled with great care.

Agents required for cleaning:

• Detergents

The cleaning has been validated using an enzymatic pH neutral detergent (Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner). An alkaline detergent might reduce the lifetime of the device.

Detergents should be used at the concentration, temperature and for the duration recommended by the manufacturer.

To remove physiological liquid inside the instrument or the spray tube, use tap water or "Aquacare" spray from Bien-Air Surgery SA.

CAUTION

Do not use detergents that are corrosive or contain chlorine, acetone or bleach, aldehydic products or alcohols.

• Lubricant:

Exclusively use "Lubrifiuid" from Bien-Air Surgery SA.

• Brush / cleaning gun

Brushes should be non-aggressive to avoid damaging the device. Preferably use nylon brushes with flexible or soft bristles. Use the cleaning gun with the appropriate nozzles for the various ducts.

If you have any additional questions regarding reprocessing procedures, instructions for reusable devices etc., please contact the regional manufacturer. Contact information can be found on the last page of these instructions.

Before use, please comply with the present section:

- Clean, lubricate and sterilize the handpiece before first use.
- Clean, lubricate and sterilize the handpiece before each further use.
- After each use, clean and sterilize the handpiece as quickly as possible.

⚠ CAUTION

Never submerge the handpiece in disinfection solutions or physiological salt water solutions (NaCl solution), as prolonged contact may cause corrosion.

9.1. Point of Use cleaning:

⚠ CAUTION

This should be performed as soon as possible after the surgical procedure.

⚠ CAUTION

Point of use cleaning must be followed by pre-cleaning.

This operation is important in order to facilitate subsequent cleaning stages (it prevents dirt from drying and sticking to the equipment).

As soon as you have finished using the handpiece, proceed as follows:

- Remove and treat the surgical instrument in accordance with the local regulations in force as regards the disposal of contaminated waste.
- Separate the handpiece from the micromotor.

For «Heavy Duty» handpieces:

- Close the locking ring.
- Rub the outside surface of the handpiece with non-woven towelettes (pre-soaked in water).

- Spray the inside of the irrigation tube in the direction it is used with “*Aquacare*” (**Fig.8**) or inject tap water with a syringe. Check that water can flow through the duct. If obstructed use a wire-pick to unblock the canal.
- Open the locking ring, remove the nose/cranio-guard and keep it separate up to the end. Rub the outside of the handpiece and the inside of the nose/cranio-guard with a non-woven towelette (pre-soaked in water).

For slim handpieces:

- Leave the locking mechanism in the open position.
- Remove the irrigation tube from the handpiece and rub it with non-woven towelettes soaked in water.
- Spray the inside of the irrigation tube with “*Aquacare*” or inject tap water with a syringe into the irrigation duct. Check that water can flow through the duct. If obstructed, use a wire-pick to unblock the canal.
- Rub the outside of the handpiece with non-woven towelettes soaked in water.

9.2. Pre-Cleaning

For «Heavy Duty» handpieces:

- Close the locking ring.
- Rinse the handpiece and nose/cranio-guard under running tap water (cold, max 20°C / 68°F) and brush the outside with nylon soft bristles for at least 30 seconds (**Fig. 9A**). Continue brushing until soil is no longer visible.
- At the end, pull and turn the locking ring to open it and brush around it with nylon soft bristles under running tap water (cold, max 20°C / 68°F) for at least 5 seconds.
- Spray (1 sec.) the inside of the irrigation tube with “*Aquacare*” in the direction it is used to remove physiological salt water.
- Perform manual or automated cleaning.

For slim handpieces:

- Close the locking ring.
- Rinse the handpiece under running tap water (cold, max 20°C / 68°F) and brush the outside with nylon soft bristles for at least 30 seconds (**Fig. 9A**). Continue brushing until soil is no longer visible.
- At the end, pull and turn the locking ring open and brush around it with nylon soft bristles under running tap water (cold, max 20°C / 68°F) for at least 5 seconds.
- Rinse the irrigation tube under running tap water (cold, max 20°C / 68°F) and brush the outside with nylon soft bristles for at least 15 seconds. Continue brushing until soil is no longer visible.
- Spray (1 sec.) the inside of the irrigation tube with “*Aquacare*” in the direction it is used to remove physiological salt water solution from the duct.
- Perform manual or automated cleaning.

9.3. Manual cleaning

For Irrigation tubes:

Automatic cleaning is not permitted.

- Soak the irrigation tube in an ultrasonic detergent solution bath for at least 3 minutes.
- Spray the inside of the irrigation tube with “*Aquacare*” in the direction it is used to remove the detergent solution.
- Thoroughly rinse under running tap water (use critical water for neurosurgical application *refer to AAMI TIR34:2014 “Water for the reprocessing of medical devices.”*) for 10 seconds. Sprays dry the irrigation duct with an airgun.

For handpieces and cranio-guards:

⚠ CAUTION

Do not place the handpiece or the cranio-guard in an ultrasonic bath

- Soak the handpiece and nose/cranio-guard in a detergent solution for at least 5 minutes.

Shake it, close and open the locking ring at least 3 times.

- Thoroughly brush the outside with nylon soft bristles for at least 30 seconds (**Fig. 9A**). Continue brushing until soil is no longer visible.
- Then, turn and release the closed locking ring and brush around it with nylon soft bristles for at least 5 seconds.
- Using a 3mm nylon brush diameter, thoroughly brush back and forth inside the handpiece with at least 3 full movements (**Fig. 9B**).
- For the cranio-guard, brush the front section thoroughly with a 3mm nylon brush, using at least three full back-and-forth movements.

Remove the handpiece and the nose/cranio-guard from the bath and thoroughly rinse vertically under running tap water (use critical water for neurosurgical application *refer to AAMI TIR34:2014 "Water for the reprocessing of medical devices."*) from both ends for at least 30 seconds moving the locking ring several times (**Fig. 10A**).

- At the end, leave it open. Make sure water goes through the inner canal. Using a syringe to inject water (critical water for neurosurgical application *refer to AAMI TIR34:2014 "Water for the reprocessing of medical devices."*) in the irrigation canal of the heavy duty handpiece.
- Dry it by wiping it with a clean and dry non-woven towelette (**Fig. 10B**). Sprays dry the inside of the handpiece from each end and the irrigation duct with an airgun.
- If the handpiece is not immediately sterilized, perform dynamic drying under ventilation, at about ~100°C (121°F), for a minimum of 25 minutes.

Or

9.4. Automatic cleaning

- Brush around the locking ring with nylon soft bristles for at least 10 seconds under running tap water. Continue brushing until soil is no longer visible on the outside. Then, close and open it at least 2 times always under running tap water. At the end, leave it open.
- Rinse the inside of the handpiece under running tap water on the coupling side for at least 15 seconds (**Fig. 10A**).
- Using a 3mm nylon brush, brush the inside back and forth from the nose side with at least 3 full movements.
- For the cranio-guard, brush the front section thoroughly with a 3mm nylon brush, using at least three full back-and-forth movements.

Place the disassembled instruments in the appropriate washer/disinfector basket and treat via a standard instrument washer/disinfector cycle (place the handpiece upright, with the tip downward).

Only use a validated washer/disinfector.

Stages of the automatic washer/disinfector cycle:

Pre-wash:

Cold tap water (<45°C / 113°F) for minimum 2 minutes.

Washing:

Hot tap water, 50°C to 60°C (122°F to 140°F) with enzymatic pH neutral detergent (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner), minimum 5 minutes.

Neutralization:

Cold tap water (<45°C / 113°F) for minimum 2 minutes.

Rinsing:

Cold critical water acc to AAMI TIR34 (<45°C / 113°F) for minimum 2 minutes.

Thermal disinfection by rinsing:

Hot critical water acc to AAMI TIR34 (90°C / 194°F) for minimum 5 minutes. Operator is responsible for the implemented value A0 according A0 concept described in EN ISO 15883 (For example, A0 600 90°C (194°F) /1min.).

Ventilated dynamic drying:

70°C (158°F), for minimum 22 minutes.

Comments:

- Comply with the washer/disinfector loading instructions provided by the manufacturer.
- Make sure that all instruments have been correctly attached to the baskets.
- Make sure that the instruments do not touch one another and that internal canals are properly rinsed.
- Remove the instruments from the washer or disinfector immediately after the machine stops and quickly lubricate and sterilize to avoid corrosion.

9.5. Inspection, lubrication and testing

Carefully inspect each part to make sure that all visible contamination has been removed. Check in particular that the ducts are clear. Where there is contamination, repeat the full cleaning process. Check for the O-ring joint on the nose/cranio-guard (**Fig. 11**) of the «Heavy Duty» handpiece.

After each cleaning operation and before each sterilization, lubricate the instrument with "Lubrifiuid" from Bien-Air Surgery SA.

To absorb any excess lubricant, apply a cloth over the instrument's apertures.

For slim handpieces, insert the specific PM2 adapter REF 1306162 of the "Lubrifiuid" spray into the back of the handpiece and squirt for about 1 second.

For «Heavy Duty» handpieces, insert the specific PM2 adapter REF 1306162 of the "Lubrifiuid" spray into the back of the handpiece and squirt for about 1 second (**Fig. 12A**) then insert it on the

front of the handpiece and squirt for about 1 second (**Fig. 12B**).

For the cranio-guard, insert the cylindrical end of the "Lubrifiuid" spray into the back of the cranio-guard and squirt for about 1 second (**Fig. 13**).

Then leave the handpiece and the nose/cranio-guard to rest upright for at least 5 minutes on a cloth.

Check the action of moving parts: Insert a rotary instrument onto the instrument and clamp the tool with the ring (**Fig. 4B**). Check that the rotary instrument remains in place when pulled.

Holding the rotary instrument between your thumb and index finger, spin the handpiece. The handpiece shall spin freely (more than 3 turns). If it does not, there is danger of burning the patient and you should send the handpiece to your retailer or Bien-Air Surgery to be repaired.

Remove the rotary instrument and leave the locking mechanism in the open position (**Fig. 3**).

Packing for sterilization:

⚠ CAUTION

Do not reassemble the handpiece and nose / cranio-guard prior to sterilization.

Separate packing: Immediately insert the handpiece in individual wrapping, such as a paper/plastic pouch or sterilization wrap for steam sterilization.

Or

Pack in stiff boxes and trays with defined, pre-configured lids and apertures and wrap the stiff boxes or tray.

Bien-Air Surgery SA can provide a sterilization tray designed for this type of item.

In the USA, FDA approved sterilization wraps, pouches or containers must be used.

9.6. Sterilization

Sterilization by steam / moist heat is recommended. The following sterilization parameters, using a Pre-Vac cycle, are recommended by Bien-Air Surgery SA to provide a sterility assurance level (SAL) of 10⁻⁶.

Only legally marketed, FDA cleared sterilizers, sterilization wrap/pouches, biological indicators etc. should be used by the end-user for packaging terminally sterilized devices.

Temperature	132°C (269°F)	134°C ² (273°F)	134°C ² (273°F)	135°C (275°F)
Time	4 min.	3 min.	18 min. ¹	3 min.
Minimum drying time	40 min. ³			

For USA only: Use sterilization cycles consistent with the cycle specifications in ANSI/AAMI ST79 "Comprehensive guide to steam sterilization and sterility assurance in health care facilities".

¹ Parameters recommended by the World Health Organization for treating instruments in the event of contamination by Non-Conventional Transmissible Agents (NCTA). Parameters recommended by Bien-Air Surgery SA.

² Not for users in US healthcare facilities.

³ Refer to sterilizer manufacturer recommendations for drying times per load configuration.

Since no reprocessing methods have been validated for removing transmissible spongiform encephalopathy (TSE) agents from medical devices, this device should not be used for patients with known or suspected TSE agent disease, including CJD and vCJD. Bien-Air Surgery recommends incinerating devices that have come into direct contact with patients suspected of or confirmed as having TSE/CJD.

The sterilizer manufacturer instructions concerning operation and load configuration should be complied with explicitly.

Do not exceed a temperature of 138°C (280°F).

Before reuse, instruments must be allowed to cool to ambient conditions, without forced cooling.

⚠ CAUTION

Never rinse instruments in cold water to cool them.

9.7. Storage

Bien-Air Surgery SA strongly advises storing only sterilized devices so as to reduce the risk of corrosion.

Ambient conditions for storage after sterilization

- Store the equipment in a clean, dry place at an ambient temperature (10–30°C (50–86°F), 20–80% humidity).
- Do not expose the equipment to direct sunlight.
- Do not expose the equipment to permanent X-ray irradiation.
- Do not store the equipment in places that could be subject to liquid splashes.
- Do not store the equipment in the following ambient conditions: - Dust - Saline or sulfurous atmosphere
- Do not store the equipment in a location where there is a risk of the release of flammable gases.

Shelf life of sterilized instruments

The shelf life of stored sterilized instruments depends on the type of packaging used and the storage conditions (refer to the DIN 58953 standard, section 9, or existing local regulations).

10. MAINTENANCE

No handpiece components are to be changed by the user. Never disassemble the handpiece.

For all servicing and repairs, we recommend that you contact your dealer or Bien-Air Surgery SA directly.

Bien-Air Surgery SA invites users to have their dynamic instruments checked or serviced once to twice a year depending on usage frequency.

Hygiene

For the safety of the repair center's personnel, the instrument should be fully cleaned and sterilized before being returned for repair. If that proves impossible, for example because disinfection or sterilization would make the instrument completely unusable, clean the instrument as carefully as possible and mark it accordingly to indicate that it has not been decontaminated.

11. MALFUNCTIONS AND ERRORS

Use the table below to solve any problems encountered. If the problem cannot be solved, stop using the product and contact a repair center approved by Bien-Air Surgery SA.

Problem	Solution
There is no transmission of movement	<p>Check that the motor operates by itself. Check that there is a rotary instrument in the locking mechanism and that it is in the "clamp" position. Make sure that the cable is properly connected to the motor and the control unit.</p> <p>Make sure the control unit is on. Check the condition of the motor power cable (broken contact or severed cable).</p>

Problem	Solution
Unable to insert the rotary instrument into the instrument	<p>Check that the chuck mechanism is in the open position (Fig. 3).</p> <p>Check that no parts are jammed inside the clamp and that the latter is clean. Check that the stem of the rotary instrument is in sound condition and of the correct diameter.</p>
Unable to close the chuck mechanism.	Manually turn the rotary instrument and close the chuck mechanism.
The rotary instrument does not rotate freely	<p>Check that the locking mechanism ring is fully inserted (Fig. 4B).</p> <p>Check that the part is clean and lubricated with "Lubrifiuid".</p>
Unable to insert the fixed cranio-guard onto the handpiece	Check that the mark on the cranio-guard is well aligned with the irrigation canal (Fig. 6C).
Excessive vibration	Stop the motor and pull on the rotary instrument to check that it is correctly in position, or change / reduce the speed of rotation.
The irrigation pump operates, but the irrigation fluid does not flow	<p>Check that the irrigation line is properly connected and that the flow rate setting is adequate.</p> <p>Check that the handpiece duct is not clogged. This should be done without a tool.</p>

12. GENERAL TERMS OF GUARANTEE

12.1. Generalities

Bien-Air Surgery SA endeavors to provide its customers with products and devices of impeccable quality, which it guarantees within the limits of the present general terms and the particular agreements signed, against any operating fault, material or manufacturing defect.

The guarantee period is 12 months from the date of invoice.

In general, the guarantee does not exempt the customer from the obligation of obtaining information from Bien-Air Surgery SA in case of doubt and in particular when the product is used in conditions that are not explicitly originally provided for.

The buyer is obliged to check the goods received within 8 days of receipt. If the goods are not checked within the period mentioned above, the customer shall be deemed to have accepted the goods, barring hidden defects. The defect notice must be received in writing by Bien-Air Surgery SA within the aforementioned period and must contain the customer's name, the date of purchase, and the product reference and serial number.

In the event of claims, Bien-Air Surgery SA or its authorized representative shall perform product repair or replacement free of charge, after analyzing the justification for the claim.

All other claims of any kind, and in particular claims for damages, are excluded.

Bien-Air Surgery SA shall not be held responsible for damage or injury and the consequences thereof, resulting in particular from:

- excessive wear;
- inappropriate use;
- failure to comply with operating instructions, assembly instructions or maintenance instructions;

- exceptional environmental, chemical, electrical or electrolytic influences;
- faulty air or water seals or electrical connections.

In any case the guarantee becomes null and void in the event of inappropriate servicing, use of non-recommended parts, accessories or consumables, or modifications to the product carried out by third parties not authorized by Bien-Air Surgery SA.

In case of dispute as to whether or not the defect exists, it shall be incumbent on the customer to prove the existence of the defect.

Guarantee claims shall be taken into consideration only upon presentation, with the product, of a copy of the invoice or delivery slip on which the date of purchase, product reference and serial number should appear clearly.

12.2. Governing law

Swiss domestic law ("Code des obligations") shall be applicable in addition to the general terms and particular agreements between the customer and Bien-Air Surgery SA.

Jurisdiction CH-2340 Le Noirmont, Switzerland.

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