REF 1600685

Rx Only

REF 2100392-0002/2022.10

**Instructions for use**  
**OSSEOSTAP**

Please read carefully and keep  
this copy for further use

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
## 1. SYMBOLS USED

	<b>CE Marking</b> with the number of the notified body.
Rx Only	<b>Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.</b>
	<b>CAUTION!</b> Refer to the accompanying documents
	<b>Consult the accompanying documents</b>
	<b>Materials to be recycled</b> The disposal and/or recycling of materials must be performed in accordance with the directives and the legislation in force.
	<b>Sterilizable in autoclave up to the specified temperature</b>
	<b>Cleaning in machine authorized</b>
	<b>Electrical or electronic materials to be recycled.</b>
	<b>Manufacturer</b>

## 2. INTENDED USE

The **OSSEOSTAP** is an active medical device that converts electrical energy into mechanical rotation to a bur during ear surgery intended for the light drilling of bones.

## 3. INDICATION FOR USE

The **OSSEOSTAP** system is intended to be used by professional users. It has been designed for the light drilling of bones of the middle ear as part of surgical operations such as stapedotomy or ossiculoplasty.

## 4. CONTRAINDICATION

None currently known.

## 5. PRECAUTIONS - WARNINGS

### 5.1. Warnings, Precautions for use

For additional information, please contact Bien-Air Surgery SA at the address indicated on the first page of this document.

**⚠ CAUTION**  
**Avoid excessive pressure on the tool. The use of a tool with excessive pressure may cause an inordinate amount of heat buildup resulting in a thermal injury to tissue. See tools instruction for use for further information.**

**⚠ CAUTION**  
**Do not use this device in the presence of a flammable anesthetics. Avoid potential ignition or explosion of gases.**

The device and its accessories should be used only by duly trained and competent medical personnel, in particular in compliance with the legal provisions in force regarding occupational safety, health and accident prevention measures, and the present user manual. According to these measures, the user has the following obligations:

- To only use devices that are in perfect working condition. In the event of irregular operation, excessive vibrations, abnormal overheating or other signs suggesting malfunctioning of the device, work must be suspended immediately. In this case, contact a repair center approved by the legal manufacturer.
- Make sure that the device is used only for the purpose for which it is intended. Protect yourself, patients and third parties from all danger and avoid contamination by the product.

The device and its accessories are designed solely for medical treatment. Any use that does not comply with its intended use is unauthorized and may prove dangerous. This medical device complies with European legal provisions in force.

Excessive noise from the bur when drilling close to the cochlea or ossicular chain may cause hearing damage.

Excessive force may lead to excess vibration, bur bending or breakage. Should a bur fracture in use, extreme care must be exercised to ensure that all fragments of the bur are retrieved and removed from the patient. Unremoved bur fragments may cause tissue damage to the patient.

Do not modify this equipment without authorization from the manufacturer. If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

Use only original maintenance products, accessories and/or spare parts approved by the legal manufacturer. The use of other products, accessories or parts could void the guarantee and/or endanger the patient or the operator.

## 5.2. Environmental protection and indications for device disposal

This equipment must be recycled. The disposal and/or recycling of materials must be performed in accordance with the directives and the legislation in force. Electrical and electronic equipment may contain dangerous substances which constitute health and environmental hazards.

The user can return the device to their distributor or directly call a firm certified to treat and recover this type of equipment (European directive 2012/19/EU).

## 6. TECHNICAL DATA

### Environmental conditions:

<b>Work</b>	Temperature:	+10°C to +30°C (+50°F to +86°F)
	Relative humidity:	20% to 80%, including condensation
	Atmospheric pressure:	700 hPa to 1060 hPa
<b>Transport</b>	Temperature:	-25°C to +70°C (-13°F to +158°F)
	Relative humidity:	10% to 100%, including condensation
	Atmospheric pressure:	600 hPa to 1060 hPa
<b>Storage</b>	Temperature:	0°C to +40°C (+32°F to +104°F)
	Relative humidity:	10% to 90%, including condensation
	Atmospheric pressure:	600 hPa to 1060 hPa

**OSSEOSTAP reference:**  
REF. 1600685

### Type:

Direct current, brushless and sensorless, curved nose.

### Power / Motor torque:

	<b>OSSEODUO</b>	ORIGO
	<b>OSSEOSTAP control unit</b>	
<b>OSSEOSTAP</b>	Max. 4W @ 3 Nmm	5 mNm

### Rotation speed:

12'000 rpm.

### Period of operation:

To avoid overheating leading to burns (temperatures between 41°C (106°F) and 48°C (118°F) on the outer surface of the **OSSEOSTAP**), the following rules must be respected:

1. Limit the maximum speed of rotation according to the bur diameter as indicated on the packaging and/or in the present operating instructions.
2. Duty cycle:  
For normal operating temperature room (typically 20°C) the micromotor are rated as below:

	<b>OSSEOSTAP</b>
<b>OSSEODUO</b>	Continuous operation
<b>OSSEOSTAP control unit</b>	
ORIGO	

For an operating room at 30°C, the micromotor are rated as below:

	<b>OSSEOSTAP</b>
<b>OSSEODUO</b>	3 min. ON, 20 min. OFF
<b>OSSEOSTAP control unit</b>	1 min. ON, 3 min. OFF
ORIGO	Continuous operation

### Noise level:

As per ISO 11498: < 60 dBA at 45 cm.

### Size / Weight:

Ø 15.5 mm, total length 180 mm / 60 g without cable and without tool.

### Bur guide tube:

Ø 2.1, length 64.5 mm, angle 15°, reinforced for optimal stability.

### Tools:

- Tungsten carbide burs Ø 0.5mm to Ø 2.3mm
- Diamond burs Ø 0.6mm to Ø 2.3mm
- Stainless steel perforator Ø 0.35mm

With color ring for diameter coding according to the table REF 1304358 supplied with the **OSSEOSTAP**.

All tools are shipped sterile packed (1 pc per box). Reusable, sterilization in autoclave.

### Motor cable:

Push-pull type connector, length 295 cm ± 5 cm.

### Compatibility:

**OSSEOSTAP** is compatible with **ORIGO**, **OSSEODUO** and **OSSEOSTAP** systems.

For further information concerning microsurgery instruments and their accessories, please contact your Bien-Air Surgery SA local distributor or consult our homepage [www.bienair.com](http://www.bienair.com).

## 7. USE

The **OSSEOSTAP** is designed to drive various types of burs. It is connected to the control unit by a specific cable.

### 7.1. Type of tools

Refer to the tools documentation. The choice of the correct instrument, especially the correct drilling diameter for stapedotomy prosthesis, is the responsibility of the surgeon.

## 7.2. Tool connection and removal

### CAUTION

Never insert an instrument on a rotating OSSEOSTAP.

1. Gently introduce the tool into the OSSEOSTAP up to the stop.
2. Completely press button A (Fig. 1) with thumb.
3. Keep button A (Fig. 1) pressed, slowly turn tool and continue to insert until white bushing C (Fig. 1) disappears.
4. Release button A (Fig. 1) and pull the bur slightly to check that it is correctly attached to the OSSEOSTAP.
5. To remove the tool, press button A (Fig. 1) and gently pull on the tool.

### CAUTION

Do not bend the bur guide tube B (Fig. 1) on the OSSEOSTAP, otherwise it can generate excessive heat during the operation.

### CAUTION

Inadequate bending of the bur guide tube can generate excessive heat during OSSEOSTAP operation.

### CAUTION

Only Bien-Air Surgery SA burs are to be used.

## 8. CLEANING / INSPECTION / STERILIZATION

### CAUTION

Do not place the OSSEOSTAP in an ultrasonic bath.

### CAUTION

Never rinse the instruments in cold water to cool them.

### CAUTION

Never submerge the OSSEOSTAP in disinfection solutions or physiological salt water solutions (NaCl solution), as prolonged contact may cause corrosion.

### CAUTION

Chemical disinfection of the motors is not recommended due to possible negative effects on the lifetime of the devices and possible residues of disinfectants.

### CAUTION

OSSEOSTAP is delivered "non-sterile".

Before it is used, please comply with the present section:

- Clean, lubricate and sterilize the OSSEOSTAP before first use.
- Clean, lubricate and sterilize the OSSEOSTAP before each further use.
- After each use, clean and sterilize the OSSEOSTAP as quickly as possible.

### Precautions for use:

Hospital procedures must be followed.

Hospital personnel working with contaminated or potentially contaminated medical instruments must comply with universal precautions.

Pointed and sharp instruments should be handled with great care.

### Agents required for cleaning:

#### • Detergents

The cleaning of the OSSEOSTAP has been validated using an enzymatic pH neutral detergent (Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner). An alkaline detergent might reduce the lifetime of the device.

Detergents should be used at the concentration, temperature and for the duration recommended by the detergent's manufacturer.

To remove physiological liquid inside the instrument or the spray tube, use "Aquacare" from Bien-Air Surgery SA.

### CAUTION

Do not use detergents that are corrosive or contain chlorine, acetone or bleach, aldehydic products or alcohols.

#### • Lubricant:

Only use "Lubrifiuid" REF 1600064 from Bien-Air Surgery SA.

#### • Brush / cleaning gun

The brushes should be non-aggressive to avoid damaging the device. Preferably use nylon brushes with flexible or soft bristles. Use the cleaning gun with the appropriate nozzles for the various ducts.

### Cleaning precautions:

Complete cleaning can lower the initial microbial load, eliminate organic matter and prevent biofilm formation. This stage is essential and influences the quality of the disinfection procedure as a whole. Cleaning combines the physico-chemical action of the product and the mechanical action of brushing and rinsing.

If you have any additional questions regarding reprocessing procedures, instructions for reusable devices etc., please contact Bien-Air Surgery SA.

### 8.1. Point of Use cleaning:

#### CAUTION

Initial cleaning should be performed at the point of use and as soon as possible after the completion of the surgical procedure.

#### CAUTION

Point of use cleaning must be followed by manual or automatic cleaning.

This operation is important in order to facilitate subsequent cleaning stages (it prevents dirt from drying and sticking to the equipment).

As soon as you have finished using the **OSSEOSTAP**, proceed as follows:

1. Disconnect the motor cable from the control unit.
2. Separate the instrument and treat it in accordance with its instructions.
3. Separate the nose from the motor: Unscrew the front part (nose) from the **OSSEOSTAP** in order to clean it separately from the motor and cable. Do not use a tool which could damage the nose.
4. Rub the outside surface of the **OSSEOSTAP** with non-woven towelettes (pre-soaked in water).
5. Spray the inside of the nose in the direction it is used with "AQUACARE" or inject tap water with a syringe.
6. Ensure that the instruments do not dry before manual or automatic cleaning by wrapping all parts in non-woven towelettes (pre-soaked in water).
7. Manual or automatic cleaning must be done max 4h after the point of use.

## 8.2. Manual cleaning

Nose:

1. Soak the nose in an enzymatic pH neutral detergent (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner) solution for at least 5 minutes.
2. While immersed, push the button for at least 3 times.
3. With a nylon soft bristles brush, brush the outside surfaces during at least 30 seconds or until soil is no longer visible.
4. With 2mm nylon brush diameter brush the inside of the nose during at least 30 seconds or until soil is no longer visible.
5. Rinse it thoroughly for at least 30 seconds under tap water with the distal part downwards.

Drying:

6. Dry it by wiping it with a clean and dry non-woven towelette. Blow dry the inside from the proximal end with pressurized air.
7. If the nose is not immediately sterilized, perform dynamic drying under ventilation, at about ~90°C (194°F), for a minimum of 25 minutes.

### OSSEOSTAP:

Hold the motor with the distal part directed downwards during the whole cleaning procedure to avoid the ingress of liquids into the motor.

Thoroughly rinse the motor and its cable under running tap water as follows:

1. Rinse the motor including the first 10cm of the cable under running tap water (cold) and brush with nylon soft bristles for at least 30 seconds on the outside surfaces. Continue brushing until soil is no longer visible.
2. Rinse the entire cable under running tap water (cold) and wipe with a wet non-woven towelette for at least 30 seconds until soil is no longer visible.

Thoroughly wash the motor and its cable with an enzymatic neutral pH detergent (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner) solution as follows:

3. Rub the entire cable with detergent impregnated non-woven wipes, moving back and forth for at least 30 seconds.
4. Wash the motor by pouring detergent solution over the device including the first 10cm of the cable, making sure to cover the complete surface.
5. Then, with a nylon soft bristles brush soaked in detergent solution, brush the outside surfaces including the first 10cm of the cable (especially in the corner, edges and openings) for at least 30 seconds. Continue brushing until soil is no longer visible.

Thoroughly rinse the motor and its cable under running tap water as follows:

6. Rinse the whole cable for at least 30 seconds.
7. Rinse the motor for at least 30 seconds.

Drying:

8. Dry the motor and the cable by wiping it with a clean and dry non-woven towelette.
9. If the micromotor and the cable are not immediately sterilized, perform dynamic drying under ventilation, at about ~90°C (194°F), for a minimum of 25 minutes.

**Or**

### 8.3. Automatic cleaning

Pre cleaning:

1. Hold the **OSSEOSTAP** and the nose with the front end directed downwards. Rinse the nose and the **OSSEOSTAP** including the first 10cm of the cable under tap water (cold, max 20°C / 68°F).
2. Clean the outside with water and a nylon brush for at least 30 seconds or until soil is no longer visible.
3. Pre-cleaning must be directly followed by final automatic cleaning.

Automatic cleaning:

Place both parts of the **OSSEOSTAP** with its cable in the appropriate washer/disinfector basket and treat via a standard instrument washer/disinfector cycle (place the nose and the motor upright, with the tip and shaft downward).

Only use a validated washer / disinfector (ISO 15883).

This product has been validated with an automated cycle consisting of the following stages:

**Pre-wash:**

Cold tap water (<45°C / 113°F) for minimum 2 minutes.

**Washing:**

Hot tap water, 50°C to 60°C (122°F to 140°F) with enzymatic pH neutral detergent (like Steris Prolystica® 2X Concentrate Enzymatic Presoak and Cleaner), minimum 5 minutes.

**Neutralization:**

Cold tap water (<45°C / 113°F) for minimum 2 minutes .

**Rinsing:**

Cold critical water acc to AAMI TIR34 (<45°C / 113°F) for minimum 2 minutes.

**Thermal disinfection by rinsing:**

Hot critical water acc to AAMI TIR34 (90°C / 194°F) for minimum 5 minutes.

Operator is responsible for the implemented value A0 according A0 concept described in EN ISO 15883 (For example, A0 600 90°C (194°F) /1min.).

**Ventilated dynamic drying:**

70°C (158°F), for minimum 22 minutes.

**Comments:**

- Comply with the washer/disinfecter loading instructions provided by the equipment's manufacturer.
- Make sure that all instruments have been correctly attached to the baskets.
- Wind up the cable correctly and without any tension (minimum inner diameter 8 cm).
- Make sure that the instruments do not touch one another and that internal canals are properly rinsed.

- Remove the instruments from the washer or disinfecter immediately after the machine stops and quickly lubricate and sterilize to avoid corrosion.

**8.4. Inspection, lubrication and testing**

Carefully inspect each part to make sure that all visible contamination and moisture has been removed. Check in particular that the ducts are clear. Where there is contamination, repeat the full cleaning process. Where there is moisture, use the airgun to dry it.

After each cleaning operation and before each sterilization, lubricate the instrument with "Lubrifiuid" from Bien-Air Surgery SA as follows:

- To absorb any excess lubricant, apply a cloth over the OSSEOSTAP's apertures.
- Spray into the drilling at the distal part of the motor shaft with "Lubrifiuid" spray for approx. 0.1 seconds.
- Spray into the nose from the proximal part with "Lubrifiuid" spray for approx. 0.1 seconds.

**⚠ CAUTION**

**Inadequate lubrication can generate excessive heat of the bur guide tube during OSSEOSTAP operation.**

- Gently reassemble the two parts of the OSSEOSTAP. Firmly tighten both parts without a tool.
- Press several times on the rotary instrument insertion / removal button in order to distribute the lubricant.
- Press several times on the rotary instrument insertion/removal button. It must always return to its initial position.

Check the device including its cable and connectors for any visible signs of deterioration or damage.

**Packing for sterilization:**

Wrap the instrument with its cable immediately after cleaning in an individual packing such as a paper/plastic pouch or sterilization wrap, for steam sterilization.

Make sure that the cable does not touch the sides of the sterilizer.

In the USA, FDA approved sterilization wraps or pouches must be used.

**8.5. Sterilization**

Sterilization by steam/moist heat is the recommended method for the OSSEOSTAP micromotor and its cable. The following table of sterilization parameters, using a B-type cycle with Pre-Vac, is recommended by Bien-Air Surgery SA to provide a sterility assurance level (SAL) of 10<sup>-6</sup>.

Only legally marketed, FDA cleared sterilizers, sterilization wrap/pouches, biological indicators etc. should be used by the end-user for packaging terminally sterilized devices.

Temperature	132°C (270°F)	134°C <sup>2</sup> (273°F)	134°C <sup>2</sup> (273°F)	135°C (275°F)
Time	4 Min.	3 Min.	18 Min. <sup>1</sup>	3 Min.
Minimum drying time	30 Min <sup>3</sup>			

<sup>1</sup> Parameters recommended by the World Health Organization for treating instruments in the event of contamination by Non-Conventional Transmissible Agents (NCTA).

<sup>2</sup> Not for users in US healthcare facilities

<sup>3</sup> Refer to sterilizer manufacturer recommendations for drying times per load configuration.

After sterilization, let the device cool down to room temperature without forced cooling.

For USA only: Use sterilization cycles consistent with the cycle specifications in ANSI/AAMI ST79 "Comprehensive guide to steam sterilization and sterility assurance in health care facilities".

Since no reprocessing methods have been validated for removing transmissible spongiform encephalopathy (TSE) agents from medical devices, this device should not be used for patients with known or suspected TSE agent disease, including CJD and vCJD. The legal manufacturer recommends incinerating devices that have come into direct contact with patients suspected of or confirmed as having TSE/CJD.

The sterilizer manufacturer instructions concerning operation and load configuration should be complied with explicitly.

**⚠ CAUTION**

**Do not exceed a temperature of 138°C (280°F).**

**⚠ CAUTION**

**Never rinse instruments in cold water to cool them.**

**8.6. Storage**

The legal manufacturer strongly advises storing only sterilized devices to reduce the risk of corrosion.

**Ambient conditions for storage after sterilization:**

- Store the equipment in a clean, dry place at an ambient temperature (See Chap. 6 for further information).
- Do not expose the equipment to direct sunlight.
- Do not expose the equipment to permanent X-ray irradiation.
- Do not store the equipment in places that could be subject to liquid splashes.
- Do not store the equipment in the following ambient conditions: Dust, Saline or sulfurous atmosphere

- Do not store the equipment in a location where there is a risk of the release of flammable gases.

**Shelf life of sterilized instruments:**

The shelf life of stored sterilized instruments depends on the type of packaging used and the storage conditions (refer to the DIN 58953 standard, section 9, or existing local regulations).

**8.7. MAINTENANCE**

No component of the **OSSEOSTAP** may be changed by the user.

The cable of the **OSSEOSTAP** is assembled to be water sealed, hence it cannot be disassembled without a specific tool. If it is very dirty on the back of the **OSSEOSTAP**, only the cable cap can be removed.

**⚠ CAUTION**

**Never disassemble the OSSEOSTAP or its cable (except the nose or the cable cap for cleaning as per instructions above).**

For all servicing and repairs, we recommend that you contact your dealer or Bien-Air Surgery SA directly.

The legal manufacturer invites users to have their dynamic instruments checked or serviced once to twice a year depending on usage frequency.

**Hygiene**

For the safety of the repair center's personnel, the instrument should be fully cleaned and sterilized before being returned for repair. If that proves impossible, for example because disinfection or sterilization would make the instrument completely unusable, clean the instrument as carefully as possible and mark it accordingly to indicate that it has not been decontaminated.

**9. MALFUNCTIONS AND ERRORS**

Use the table below to solve any problems encountered. If the problem cannot be solved, stop using the product and contact a repair center approved by the legal manufacturer.

Problem	Solution
<p>The <b>OSSEOSTAP</b> is not working.</p> <p>Impossible to properly insert a tool into the <b>OSSEOSTAP</b></p>	<p>Make sure that the cable is properly connected to the control unit.</p> <p>Make sure the control unit is turned on.</p> <p>Check the condition of the motor cable and the <b>OSSEOSTAP</b> (broken contact or severed cable).</p> <p>Check the condition of the tool.</p> <p>Without a tool, check the freedom of movement of the Bur insertion / removal button A (<b>Fig. 1</b>).</p> <p>To do so, push it up to the thrust stop and release. It should resume to its initial position. Check that there is no part or other body retained inside the bur guide tube B (<b>Fig. 1</b>),</p> <p>To do so, unscrew the nose from the <b>OSSEOSTAP</b>. Remove any obstructions from the guide tube</p>

Problem	Solution
Bur insertion / removal button A (Fig. 1) is not free to move.	Unscrew the nose from the OSSEOSTAP and check cleanliness of the button.
The tool is inserted, but it does not hold in place	Lubricate with "Lubrifluid" before sterilization. Check that the tool is correctly inserted: White bushing C (Fig. 1) on tool must completely disappear. If the problem persists, contact your accredited Bien-Air Surgery SA distributor for repair.

## 10. GENERAL TERMS OF GUARANTEE

### 10.1. General Information

Bien-Air Surgery SA endeavors to provide its customers with products and devices of impeccable quality which are guaranteed, within the limits of the present general terms and any specific agreements signed, against any operating fault, material or manufacturing defect.

The guarantee period is 12 months from the date of invoice.

In general, the guarantee does not exempt the customer from the obligation to obtain information from Bien-Air Surgery SA in case of doubt and in particular when the product is used in conditions not explicitly provided for originally.

The buyer is obliged to check the goods received within 8 days following their receipt. If the goods are not checked within the aforementioned period, the customer shall be deemed to have accepted the goods, barring hidden defects. The defect notice

must be received in writing by Bien-Air Surgery SA within the aforementioned period and must contain the customer's name, the date of purchase, the product reference and serial number.

In the event of claims, Bien-Air Surgery SA or its authorized representative shall perform product repair or replacement free of charge, after analyzing the justification for the claim.

All other claims of whatsoever kind, and in particular claims for damages, are excluded.

Bien-Air Surgery SA shall not be held responsible for damage or injury and the consequences thereof, resulting in particular from:

- excessive wear,
- inappropriate use,
- failure to comply with operating instructions, assembly instructions or maintenance instructions,
- exceptional environmental, chemical, electrical or electrolytic influences,
- faulty air or water seals or electrical connections.

In any case the guarantee becomes null and void in the event of inappropriate servicing, use of non-recommended parts, accessories or consumables, or modifications to the product carried out by third parties not authorized by Bien-Air Surgery SA.

In case of dispute as to whether or not the defect exists, it shall be incumbent on the customer to prove the existence of the defect.

Guarantee claims shall be taken into consideration only upon presentation, with the product, of a copy of the invoice or delivery slip on which should appear clearly the date of purchase and the product reference and serial number.

### 10.2. Governing law

Swiss domestic law (Code of obligations) shall apply in addition to the general terms and conditions and particular agreements between the customer and Bien-Air Surgery SA.

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