

PRODUCTS FOR INTERVENTIONAL SIALENDOSCOPY

Minimally invasive solutions designed for gland preservation

Access and treatment essentials

INTERVENTIONAL SIALENDOSCOPY

*The minimally invasive treatment of
obstructive salivary gland disease*

**ACCESS
PRODUCTS →**

**INTERVENTION
PRODUCTS →**

**ORDERING
INFORMATION →**

REDUCED
RECOVERY TIME



NO
INCISION*



GLAND
PRESERVATION



Sialendoscopy offers a variety of advantages to the patient, including reduced recovery time, no incision or subsequent scarring, minimized potential for nerve damage, preservation of the salivary gland, and early return to a normal diet.^{1, 2, 3}

Sialendoscopy allows stones and strictures to be managed in office with a low risk of major complications to the patient compared to open surgical procedures.¹

*A combined approach may be necessary depending on stone size, mobility, and clinical judgement.



Cook Medical continually collaborates with physicians to improve procedures. Cook has partnered with leading physicians to create a comprehensive line of salivary duct access and interventional products that are designed to maintain access throughout the procedure and treat obstructive salivary gland disorders. This line of products is designed to help reduce trauma for the patient and increase efficiency with every step.


Cook's salivary duct access tools allow me to gain and maintain access to the duct, which gives me much greater control throughout the entire procedure.

Jack Kolenda, MD

Chief, Department of Otolaryngology
Oakville-Trafalgar Memorial Hospital
Oakville, Ontario


Dr. Kolenda is the inventor of the Kolenda Salivary Access Introducer Set and received a royalty payment based on Cook's license to use this technology.

Gain and maintain access to the salivary duct throughout the entire procedure.




Wire Guide

Access the salivary duct through the papilla.



Salivary Access Dilator Set


Introduce a series of flexible dilators over the wire guide to expand the papilla and prepare the salivary duct for the introduction of procedural instruments.



Kolenda Salivary Access Introducer Set


Create an open working channel into the salivary duct by passing the sheath over the wire guide. The sheath is also designed to protect the ductal wall as you insert and remove procedural instruments.

Gain and maintain access to the salivary duct throughout the entire procedure.




NCircle® Tipless Salivary Stone Extractor

Insert the NCircle nitinol basket into the duct to remove salivary stones and fragments. The stone extractor has a tipless configuration that allows you to position the basket against the mucosal lining.




NGage® Salivary Stone Extractor

Insert the NGage nitinol basket into the duct to remove salivary stones and fragments. The stone extractor works as both a basket and a grasper.



SialoCath® Salivary Duct Catheter

Insert the catheter over a wire guide or through the Kolenda Salivary Access Introducer to irrigate the submandibular or parotid salivary ducts to flush stone fragments.



Advance® Salivary Duct Balloon Catheter Set

Once you have created a working channel in the duct, insert the balloon catheter over a wire guide and use it alongside a sialendoscope for direct visualization to dilate at the stricture site.

ORDERING INFORMATION

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Wire Guide


Order Number	Reference Part Number	Wire Guide Diameter inch	Wire Guide Length cm	Tip Configuration
G02274	C-SF-15-50	0.015	50	Straight


Salivary Access Dilator Set

Order Number	Reference Part Number	Dilator Fr	Dilator OD mm	Dilator OD inch
G24131	SDDS-18-20-HC	4.0	1.40	0.055
		5.0	1.67	0.066
		6.0	2.0	0.079
		7.0	2.33	0.092

Kolenda Salivary Access Introducer Set

Order Number	Reference Part Number	Introducer ID mm/inch	Wire Guide Diameter inch	Sheath Length cm
G24129	SDAS-6.0-18	2.21/0.087	0.018	5.5
G24132	SDAS-5.0-18	1.88/0.074	0.018	5.5





NGage® Salivary Stone Extractor

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Basket Diameter mm
G23798	SSEB-1.7-115-8	1.7	115	8

NCircle® Tipless Salivary Stone Extractor

Order Number	Reference Part Number	Catheter Fr	Catheter Length cm	Basket Diameter mm
G23797	SSEB-1.5-115-10	1.5	115	10

SialoCath® Salivary Duct Catheter

Order Number	Reference Part Number	Catheter OD mm/inch	Catheter ID mm/inch	Compatible Wire Guide Diameter inch	Compatible Catheter Length cm	Tip Configuration
G35166	SDC-2.2-20	0.74/0.029	0.53/0.021	0.018	15	Straight

Advance® Salivary Duct Balloon Catheter Set

Order Number	Reference Part Number	Wire Guide Diameter inch	Inflated Balloon Diameter mm
G34991	SDBS-14-30-1.5-2	0.014	1.5
G34992	SDBS-14-30-2.0-2	0.014	2.0
G34993	SDBS-14-30-2.5-2	0.014	2.5
G34994	SDBS-14-30-3.0-2	0.014	3.0

REFERENCES

1. Nahlieli O, Nakar LH, Nazarian Y, et al. Sialendoscopy: a new approach to salivary gland obstructive pathology. *J Am Dent Assoc.* 2006;137(10):1394-1400.
2. Carta F, Farneti P, Cantore S, et al. Sialendoscopy for salivary stones: principles, technical skills and therapeutic experience. *Acta Otorhinolaryngol Ital.* 2017;37(2):102-112.
3. Achim V, Light TJ, Andersen PE. Gland preservation in patients undergoing sialendoscopy. *Otolaryngol Head Neck Surg.* 2017;157(1):53-57.

Proudly distributed by



Endotherapeutics
IMPROVING HEALTHCARE



📞 **1300 00 ENDO (3636)**
info@endotherapeutics.com.au
www.endotherapeutics.com.au
Suite 301, 2 Banfield Road
Macquarie Park NSW 2113

AUSTRALIAN SPONSOR:

William A. Cook Australia Pty Ltd
ABN 79 005 526 723
95 Brandl Street
Eight Mile Plains QLD 4113