


A-Cysta-Cath™

Suprapubic Catheter Introducer

The A-Cysta-Cath has been designed specifically to allow for the introduction of a suprapubic catheter with minimal time and effort.



12 Fr A-Cysta-Cath

Order code: YM SCIS-001P (Box of 10)

16 Fr A-Cysta-Cath

Order code: YM SCIS-002P (Box of 10)



Endotherapeutics

Suite 1, Level 1, 51 Rawson Street
Epping NSW 2121, Australia
T: 1300 00 ENDO (3636) F: + 61 2 9869 7723
E: info@endotherapeutics.com.au
www.endotherapeutics.com.au

A-Cysta-Cath™ Plus

SupraCath suprapubic catheter with matching introducer all in one pack!

Now you can enjoy the ease of ordering a two-way suprapubic catheter and corresponding introducer and scalpel in one convenient pack.



Why choose the SupraCath?

Short soft tip reduces bladder irritation and is ideal for long-term drainage

Open tip design prevents blockages and enables over-the-wire exchange

Large side eyes allow maximum drainage

Soft, medical-grade silicone is hypoallergenic and less prone to encrustation

NEW! Depth measurements to ensure accurate placement during catheter exchange

Integrated balloon enables insertion through a smaller introducer, facilitates smooth removal by minimising cuff difficulties, **large diameter lumen** ensures optimum drainage

18 Fr A-Cysta-Cath Plus Kit

16 Fr A-Cysta-Cath + 18 Fr SupraCath + Scalpel

Order code: YM-1034-18 (Box of 10)

14 Fr A-Cysta-Cath Plus Kit

12 Fr A-Cysta-Cath + 14 Fr SupraCath + Scalpel

Order code: YM-1034-14 (Box of 10)



Endotherapeutics

Suite 1, Level 1, 51 Rawson Street
Epping NSW 2121, Australia
T: 1300 00 ENDO (3636) F: + 61 2 9869 7723
E: info@endotherapeutics.com.au
www.endotherapeutics.com.au