

# You asked for it! Introducing the 22 Fr. Gel Cap

NEW



22 Fr. Gel Caps  
BE 3295 (60 ea.)



20 Fr. Gel Caps  
BE 3290 (90 ea.)



18 Fr. Gel Caps  
BE 3195 (90 ea.)



16 Fr. Gel Caps  
BE 3190 (90 ea.)

*(Gel caps not shown actual size)*

The 22 Fr. Blom-Singer® Replacement Gel Cap for gel cap insertion is here!

The gel cap insertion method to permit ease of prosthesis insertion has been a bedrock of the Blom-Singer® voice prosthesis product line. In response to the feedback from our valued customers, we now offer gel caps in a wide variety of sizes including 16 Fr., 18 Fr., 20 Fr. and 22 Fr.

**INHEALTH®**  
TECHNOLOGIES

# Blom-Singer® Voice Prosthesis product line and recommended gel cap sizes *(intended as a guide only)*

| Prosthesis   | 16 Fr. | 18 Fr. | Prosthesis   | 20 Fr. | 22 Fr. |
|--------------|--------|--------|--------------|--------|--------|
| DB 16 XX     | ●      |        | BE 20XX      | ●      |        |
| BE 60 XX     | ●      |        | LP 20XX      | ●      |        |
| LP 16XX      | ●      |        | LP 20XX-IR   | ●      |        |
| LP 16XX-IR   | ●      |        | LP 20XX-NS   | ●      |        |
| IN 16XX-NS   | ●      |        | LP 20XX-SO   | ●      |        |
| IN 16XX-SO   | ●      |        | LP 20XX-SL   | ●      |        |
| IN 16XX-SL   | ●      |        | LP 20XX-IR   | ●      |        |
| IN 16XX-IR   | ●      |        | AS 20XX-NS   | ●      |        |
| AS 16XX-NS   | ●      |        | AD 20XX-NS   | ●      |        |
| IN 16XX-LF   | ●      | ●      | IN 20XX-LF   | ●      | ●      |
| IN 16XX-LEF  | ●      | ●      | IN 20XX-LEF  | ●      | ●      |
| IN 16XX-LETO | ●      | ●      | IN 20XX-LETO | ●      | ●      |
| IN 16XX-LEIR | ●      | ●      | IN 20XX-LEIR | ●      | ●      |
| IN 16XX-LESL | ●      | ●      | IN 20XX-LESL | ●      | ●      |
|              |        |        | DV 20XX-NS   | ●      |        |
|              |        |        | DV 20XX-LF   | ●      | ●      |

Prescription required for all gel caps.

Proudly distributed by



1300 00 ENDO

info@endotx.com.au

endotx.com.au

Suite 301, 2 Banfield Road  
Macquarie Park NSW 2113