

INTRODUCING

The Blom-Singer® Electrolarynx EL 2000

Innovation in Every Word



The Next Generation of Digital Speech Aid

We are thrilled to introduce the Blom-Singer® Electrolarynx EL2000 by InHealth Technologies®, an innovative device designed to break communication barriers and enhance speech clarity.

Features:

- Activator Button (red)
- Pitch/volume rocker signs (+) (-)
- Built-in battery
- Up to 30 hrs speaking time
- Magnetic charging USB port
- Intraoral accessory (EL2002 Oral Adapter*) *Sold Separately

Effortless Communication

Engineered for continuous engagement, the EL2000 empowers you to participate in conversations confidently. Gain more control over your communication with the intraoral accessory, ensuring comfortable and effective use.

Designed for Comfort and Simplicity

Lightweight and ergonomically designed, the EL2000 offers comfortable and reassuring daily use. Intuitive controls allow for easy adjustments, helping you stay connected and present in every moment.

Description	Product Code	CH2 Order Code	Cost
Blom-Singer® Digital ElectroLarynx	IHT EL200	2804264	\$1100
EL2002 ORAL ADAPTER	IHT EL2002	2804272	\$75



Endotherapeutics
IMPROVING HEALTHCARE

1300 00 ENDO
info@endotx.com.au
endotx.com.au
Suite 301, 2 Banfield Road
Macquarie Park NSW 2113