


INTRODUCING.... THE OTO-RHINO BALLOON EXTRACTOR

NASAL FOREIGN BODY REMOVER

We would like to offer your Emergency Department the **enclosed FREE TRIAL BALLOON EXTRACTOR** specifically designed to remove foreign bodies from the nasal or ear passage often presented by Paediatric Patients.

They come individually packaged and are latex free. The 1.5cc stopper on the syringe is specifically adapted to prevent over-inflation. The catheter comes covered in a protective sheath. Assembly is quick and easy: simply screw the two parts together and begin use.



Easy Extraction

Allows quick and easy extraction of nasal foreign bodies of almost any type

Atraumatic

Soft silicone balloon adhered to a thin catheter provides painless removal

Ready to Use

Complete single use product comes ready without assembly or the need for additional equipment

Scan to watch the Instructional Video



Low Cost

Low price compared to other extractors on the market and other traumatic methods which may require sedation

ARTG #339035
1 Box Contains 10 units
Order Code: MAR BEX-7010

Proudly distributed by



1300 00 ENDO
info@endotx.com.au
endotx.com.au
Suite 301, 2 Banfield Road
Macquarie Park NSW 2113