


Is your Laryngectomy 'Protector' really protecting you or is it just a cover?


The Buchanan range is proven to filter up to 49% more airborne particles than the competition⁽¹⁾

The filtration ability of Buchanan® and Buchanan® Lite stoma protectors were tested against two competitor stoma covers in a laboratory⁽¹⁾. The graph below shows the percentage of particles which each product was able to filter. Airborne particles equivalent to the size of mould spores, dust and some bacteria were tested (4.5 microns). These are small enough to be inhaled and reach the middle of the lungs.




4 layers

Buchanan® Lite - with dual layer microfibre core



3 layers

Buchanan® - Cotton mesh with Hydrolox® foam core



1 layer

Stoma covers 1 and 2 - 3D mesh fabric

The Buchanan® range - designed and tested to protect your airway, not just cover it

(1) Laryngectomy Protector Filter Media Tests - Report No. 17398 - Particle Technology Ltd - 24th May 2018

Scientifically proven protection, filtration and humidification for tracheostomy and laryngectomy.



Buchanan®

Buchanan Large

Dimensions:

Width: 216 mm Length: 208 mm

LA BUP 0001	Large	Pack of 10
-------------	-------	------------

Buchanan Small

Dimensions:

Width: 140 mm Length: 114 mm

LA BUP 0002	Small	Pack of 10
-------------	-------	------------



Buchanan® Lite

Buchanan® Lite Large

Dimensions:

Width: 225 mm Length: 185 mm

LA BLT 1001	Large	White	Pack of 10
LA BLT 1002	Large	Blue	Pack of 10
LA BLT 1003	Large	Beige	Pack of 10

Buchanan® Lite Small

Dimensions:

Width: 140 mm Length: 115 mm

LA BLT 2001	Small	White	Pack of 10
LA BLT 2002	Small	Blue	Pack of 10
LA BLT 2003	Small	Beige	Pack of 10

Buchanan® DeltaNex



Buchanan® DeltaNex Protectors

LA DNX 0001	Neck straps	500 mm	Pack of 2
LA DNX 0002	Neck Straps	600 mm	Pack of 2
	600 mm	Pack of 2	