

Blind biopsies can  
leave you in the dark

Blind biopsy can miss diagnoses and may add to treatment delays



50%

If cancer occupies less than 50% of the surface area of the endometrial cavity, the cancer can be missed by a blind endometrial biopsy<sup>1</sup>



18%

Blind biopsy alone could potentially miss the diagnosis of focal lesions in up to 18% of patients<sup>2</sup>




In the evaluation of abnormal uterine bleeding (AUB)

Endosee<sup>®</sup> simply gives you direct visualization at point of care to get to a diagnosis quicker



In premenopausal ovulatory women with AUB, hysteroscopy will detect an anatomical structural lesion in

65%–80% OF PATIENTS<sup>3</sup>




Endosee allows you to diagnose uterine focal pathology more accurately than endometrial biopsy alone



Sensitivity for diagnosing polyps<sup>4</sup>



Sensitivity for diagnosing myomas<sup>4</sup>



Sensitivity for diagnosing hyperplasia<sup>4</sup>



# Endosee—the point-of-care diagnostic hysteroscope

Simple, all-in-one hysteroscope allows direct visualization of entire uterine cavity

Handheld, portable, cordless system is quick to set up in any room at any time


Sterile, single-use, thin (<5 mm), flexible cannula with camera and light source provides excellent visualization

Allows in-office diagnostic hysteroscopy to be completed in minutes, either alone or during same visit as an endometrial biopsy

## Ordering information



**ES8000**  
ENDOSEE HYSTEROSCOPE  
(includes handset, docking station, AC adapter, USB cable)

**ESDX5**  A long, thin, flexible diagnostic cannula with a camera lens at the tip.  
**ENDOSEE DX CANNULA**  
(disposable, single-use, sterile diagnostic cannula, box of 5). Includes a free endometrial sampler with every cannula.

## Optional accessories

**ESSYR**  
50mL SYRINGE  
(sterile 50/case)

**ESTUB**  
STD BORE EXT TUBING 30IN  
(sterile 100/case)

To order, contact a CooperSurgical representative.

**CooperSurgical**

©2016 CooperSurgical, Inc. 82653 5/16



Endotherapeutics

References: 1. ACOG Practice Bulletin Number 128. 2. Goldstein S, Zeltser I, Horan C, et al. Am J Obstet Gynecol . 1997;177:102-108. 3. Isaacson K. Curr Opin Obstet Gynecol . 2002;14:381-385. 4. Garcia A. OBG Manage . 2013;25:44-48.

**Distributed by Endotherapeutics**

ABN 89 088 045 103

Level 1, 51 Rawson St Epping NSW 2121 Australia

T: 1300 00 ENDO F: +61 (2) 9869 7723 E: info@endotherapeutics.com.au

[www.endotherapeutics.com.au](http://www.endotherapeutics.com.au)



**endosee**<sup>®</sup>  
office hysteroscopy