

EndoSee® hysteroscopy simplifies: portability

setup

space

cost

visualization

everything



endoSee®
office hysteroscopy


Introducing EndoSee— a truly innovative way to perform in-office, diagnostic hysteroscopy in any room, at any site

An all-in-one, handheld, portable, cordless system for diagnostic hysteroscopy

Simple and quick to set up, requiring minimal staff training

System includes a low-cost, reusable handset and an affordable, single-use cannula

Comes with dock for battery charge and data transfer via USB
(rechargeable battery lasts at least 2 hours)





Ordering information

ES8000
ENDOSEE HYSTEROSCOPE
(includes handset, docking station,
AC adapter, USB cable)

ESDX5
ENDOSEE DX CANNULA
(disposable, single-use, sterile diagnostic
cannula, box of 5)

Optional accessories

ESSYR
50mL SYRINGE (sterile 50/case)

ESTUB
STD BORE EXT TUBING 30IN (sterile 100/case)

CooperSurgical

©2015 CooperSurgical, Inc. 82593 04/15



Endotherapeutics

Distributed by Endotherapeutics
ABN 89 088 045 103
Level 1, 51 Rawson St Epping NSW 2121 Australia
T: 1300 00 ENDO F: +61 (2) 9869 7723 E: info@endotherapeutics.com.au
www.endotherapeutics.com.au



endosee®
office hysteroscopy