



**NEW**

# Innovex Flex-Tip Negative Pressure Access Sheath

Designed to improve safety and efficiency during urological procedures



**Endotherapeutics**  
IMPROVING HEALTHCARE

# Innovex Flex-Tip Negative Pressure Access Sheath

The Innovex Flex-Tip Negative Pressure Access Sheath is an innovative device designed to improve safety and efficiency during urological procedures.

The Innovex Flex-Tip Negative Pressure Access Sheath features a unique negative pressure suction system that actively removes debris and fluids from the urinary tract.



- 1. Enhanced Safety:** Reduces intrarenal pressure during ureteroscopy in the kidney to minimise the risk of sepsis and injury. Reduces the need for multiple reinsertions of the ureteroscope which can lower the risk of ureteral trauma.
- 2. Increased Efficiency:** Improves visualisation in the ureter and kidney during laser lithotripsy. Enhances stone debris and fragment retrieval which reduces overall procedure times.
- 3. Improved ease of use:** Kink-resistant hydrophilic coated sheath is reinforced with a stainless steel coil, designed to provide stability and strength through tortuous ureters and reduces friction and irritation to the ureter.
- 4. The New Innovex Flex-Tip Negative Pressure Access Sheath** allows for the access sheath to reach the target calyx and stone with the flexible ureteroscope.
- 5. Flexible application:** Our unique design allows suction control when needed, giving the clinician maximum choice in the clinical application all from one sheath.

The New Innovex Flex-Tip Negative Pressure Access Sheath can reduce intrarenal pressure, improve vision, remove stone fragments and save operating time. With its unique negative pressure system, reinforced tip, and patient-friendly design, the Innovex Flex-Tip Negative Pressure Access Sheath is poised to become the preferred choice for urologists.

Code	Size
INN IVX-NS-1040	Innovex Flex-Tip Negative Pressure Access Sheath - 10/12F x 40cm
INN IVX-NS-1050	Innovex Flex-Tip Negative Pressure Access Sheath - 10/12F x 50cm
INN IVX-NS-1140	Innovex Flex-Tip Negative Pressure Access Sheath - 11/13F x 40cm
INN IVX-NS-1150	Innovex Flex-Tip Negative Pressure Access Sheath - 11/13F x 50cm
INN IVX-NS-1240	Innovex Flex-Tip Negative Pressure Access Sheath - 12/14F x 40cm
INN IVX-NS-1250	Innovex Flex-Tip Negative Pressure Access Sheath - 12/14F x 50cm

**Contact us today to learn how the Innovex Flex-Tip Negative Pressure Access Sheath can benefit your practice and patients**

Proudly distributed by



1300 00 ENDO  
 info@endotx.com.au  
 endotx.com.au

