



# MEDIPLAST<sup>®</sup> ENT PRODUCTS



**FAST EASY SAFE**



**Endotherapeutics**

Suite 1, Level 1, 51 Rawson Street, Epping NSW 2121

Phone: 1300 00 ENDO (3636) or (02) 9869 2868

Fax: (02) 9869 7723

Email: [info@endotherapeutics.com.au](mailto:info@endotherapeutics.com.au)

[www.endotherapeutics.com.au](http://www.endotherapeutics.com.au)

# MEDIPLAST® ENT PRODUCTS

The Mediplast ENT products are the leading single-use choice.


**High Quality - Safe - Easy - Fast**

Doctors all over the world prefer the Mediplast products for numerous reasons:

- **Clean and sterile**  
It is time consuming and difficult to clean reusable suction instruments. There is a trend worldwide to use single-use instruments in healthcare to improve quality control and avoid cross-contamination.
- **Cost efficient**  
The “real” price of a reusable instrument is the cost of cleaning and re-sterilising. Compared to a disposable product, this is often higher than the cost of a Mediplast disposable product.
- **Specially designed for ENT examination and surgery**  
Mediplast metal and plastic suction tips are designed in close collaboration with ENT doctors. Important features such as the angle of the suction tip, which enables the doctor to have an unimpeded view through a microscope, and the suction control handle connector that facilitates quick change of suction tips during surgery, were all the result of ensuring the needs of an ENT Specialist were met.



# Single-use ENT Suction Devices




**Fast  
Easy  
Safe**

# OTOLOGY

- Trans/Intratympanic Injection Needle
- Metal Suction Tips
- Aural Speculums
- Suction Handles
- Medibender
- Myringotomy Scalpels
- Injection/Aspiration Canulas
- Specal Canulas
- Myringotomy Sets
- Otology Sets
- Suction Sets



## Trans/Intratympanic Injection Needle




Specially designed needle for safer and easier transtympanic injections to the middle ear. The needle is 90mm long and has a specially designed sharpened tip end which is 2mm long and holds an outer diameter of only 0,4mm. The principle of “tube in tube” is used to create a natural stop against the tympanic membrane when injecting. Designed in association with Prof Måns Magnusson of the ENT clinic at Lund University Hospital.

Code	Size	Colour	Description
6066590	0,4 x 90 mm	Lilac	Luer-lock

## Metal Suction Tip



Made of medical grade stainless steel (virgin steel)

The angle of the Suction Tip is specially shaped for unobstructed view by examination and surgery through microscope



Patient friendly  
Beveled and smooth tip end

Colour coded Luer-lock hub of polypropylene that fits on all standard Luer lock and Luer slip parts





## Metal Suction Tip

- Sterile, single use
- For suction of wax, puss, blood, secretion and others
- For dosing of salves and medicine
- For ENT surgical incisions
- For dental, mouth and microsurgery



## Metal Suction Tips

Code	Size	Gauge	Colour	Description
Medi 66500040	0.4 x 70 mm	27G	Purple	Luer lock
Medi 66507040	0.4/0.7 x 77 mm	27/22G	Purple	Luer lock
Medi 66500070	0.7 x 70 mm	22G	Black	Luer lock
Medi 66500100 *	1.0 x 70 mm	19G	Brown	Luer lock
Medi 66500120	1.2 x 70 mm	18G	Pink	Luer lock
Medi 66500140 *	1.4 x 70 mm	17G	White	Luer lock
Medi 66501140 *	1.4 x 70 mm	17G	White	Round hub
Medi 66500170 *	1.7 x 70 mm	15.5G	Grey	Luer lock
Medi 66500208 *	2.0 x 70 mm	14G	Green	Luer lock
Medi 66502140 *	2.0 x 70 mm	14G	Green	Round hub
Medi 66500045	0.4 x 80 mm	27G	Purple	Luer lock
Medi 66500072	0.7 x 80 mm	22G	Black	Luer lock
Medi 66500102	1.0 x 80 mm	19G	Brown	Luer lock
Medi 66500122	1.2 x 80 mm	18G	Pink	Luer lock
Medi 66500142 *	1.4 x 80 mm	17G	White	Luer lock
Medi 66501180 *	1.4 x 80 mm	17G	White	Round hub
Medi 66500172	1.7 x 80 mm	15.5G	Grey	Luer lock
Medi 66500200 *	2.0 x 80 mm	14G	Green	Luer lock
Medi 66502200 *	2.0 x 80 mm	14G	Green	Round hub
Medi 66500250 *	2.5 x 80 mm	12G	Orange	Round hub
Medi 66500300 *	3.0 x 80 mm	11G	Blue	Round hub
Medi 66500400 *	4.0 x 80 mm	8G	Yellow	Round hub
Medi 66509100	1.0 x 100 mm	19G	Brown	Luer lock
Medi 66509122	1.2 x 100 mm	18G	Pink	Luer lock
Medi 66509140 *	1.4 x 100 mm	17G	White	Luer lock
Medi 66501100	1.4 x 100 mm	17G	White	Round hub
Medi 66509170	1.7 x 100 mm	15.5	Grey	Luer lock
Medi 66509200 *	2.0 x 100 mm	14G	Green	Luer lock
Medi 66502100 *	2.0 x 100 mm	14G	Green	Round hub
Medi 66509250 *	2.5 x 100 mm	12G	Orange	Round hub
Medi 66509300 *	3.0 x 100 mm	11G	Blue	Round hub

\* also available in big-pack

# The items highlighted are kept in stock, others if ordered may have a short wait time on ordering





## Metal Suction Tip (Stapes)

- Sterile, single use
- Tube in tube for rigidity and fine tip end lumen
- For stapedotomy and stapedectomy procedures



## Metal Suction Tips Type Olive

Code	Size	Gauge	Colour	Description
Medi 66500510	1.0 x 60 mm	19G	Brown	Luer lock
Medi 66510070	1.0 x 70 mm	19G	Brown	Luer lock
Medi 66510080	1.0 x 80 mm	19G	Brown	Luer lock
Medi 66510100	1.0 x 100 mm	19G	Brown	Luer lock
Medi 66500514	1.4 x 60 mm	17G	White	Luer lock
Medi 6514070	1.4 x 70 mm	17G	White	Luer lock
Medi 66514080	1.4 x 80 mm	17G	White	Luer lock
Medi 66514100	1.4 x 100 mm	17G </td <td>White</td> <td>Luer lock</td>	White	Luer lock
Medi 66517070	1.7 x 70 mm	15.5G	Grey	Luer lock
Medi 66517080	1.7 x 80 mm	15.5G	Grey	Luer lock
Medi 66517100	1.7 x 100 mm	15.5G	Grey	Luer lock
Medi 66500520	2.0 x 60 mm	14G	Green	Luer lock
Medi 66520070	2.0 x 70 mm	14G	Green	Luer lock
Medi 66520080	2.0 x 80 mm	14G	Green	Luer lock
Medi 66520100	2.0 x 100 mm	14G	Green	Luer lock



## Aural Speculum

- Hygienically packed
- Maximal view via oval end
- Dull surface to prevent reflexions
- Angled cut at tip end
- Wide brim for stability
- 7 sizes, fits Siegle pneumatic head

Code	Size
Medi 606625	2.50 mm
Medi 606630	3.25 mm
Medi 606640	4.00 mm
Medi 606650	5.00 mm
Medi 606660	6.00 mm
Medi 606670	7.00 mm
Medi 606680	8.00 mm



## Siegle Pneumatic Head


Code	Description
Medi 66508	Siegle pneumatic head



## Suction Tip Handle

- Sterile, single use
- Fits all metal suction tips with Luer slip hub
- With or without vacuum control
- Allows quick change of tips during surgery
- Made of white ABS plastic
- Tip bending holes

Code	Description
Medi 6650992	With suction control
Medi 6650509	Without suction control
Medi 66500990	With suction control
Medi 66500991	With suction control



## Medibender

### A bending instrument for metal suction tips

- Stainless steel, reusable
- For easy and controlled bending of metal suction tips 1.0, 1.4, 2.0, 3.0 and 4.0 mm diameters
- For the tip end to reach tight areas like atticus where a little bend is needed to reach all the way


Code	Description
Medi 66500	For 2.0, 1.4 and 1.0 mm
Medi 66501	Also for 3.0 and 4.0 mm



## Myringotomy Scalpel

3 edged cutting edge or traditional  
4 edged spear head shape  
No need for cutting, only puncture

The scalpel is shaped to an  
angle for free sight when  
performing paracentesis



Made of medical grade stainless  
steel (virgin steel)


Handle with 2 finger grip for good  
balance and stability  
Can be connected with suction tip  
for extra stability



## Myringotomy Scalpel

- Sterile, single use
- For puncture of the tympanic membrane
- Single use - Guaranteed sharp edge for every paracentesis

Code	Description
Medi 66121070	1.0 x 70 mm pink, triple cut edges
Medi 66121370	1.3 x 70 mm red, triple cut edges
Medi 66121170	1.0 x 70 mm orange, double edged



*# The items highlighted are kept in stock, others if ordered may have a short wait time on ordering*

## Injection Canulas / Aspiration Canulas



- For numerous applications
- Canulas with Quincke grinding and thin wall for increased flow
- Translucent and color-coded luer-lock to simplify the identification
- Protective caps made of polypropylene
- Single use, sterile and single packed
- Material: PP/PE/Medical Grade Stainless Steel

Code	Size	Gauge	Colour
6063560080	0,8 x 80	21	Green
6063560205	2,1 x 80	14	Green
6062640440 (dental)	0,4 x 40	27	Grey



## Special Canulas with and without Mandrin

- Hypodermic needles with specially bevelled points
- For numerous applications.
- Single use, sterile and single packed
- With metal hub and luer-lock
- Protective caps made of polypropylene
- Material: PP/PE/Medical Grade Stainless Steel



Code	Size	Gauge	Colour
<b>Canulas, special with Quincke grinding</b>			
6062590715	0,7 x 150	22	14
6062590720	0,7 x 200	22	14
6062590811	0,8 x 110	21	14
6062600093	0,9 x 110	20	14
6062620091	0,9 x 120	20	14
6062630090	0,9 x 150	20	14
6062630192	0,9 x 200	20	14
6062630195	0,9 x 250	20	14
6062610126	1,2 x 115	18	14
6062600135	1,3 x 110	18	20
6062600184	1,8 x 110	15	20
6062610209	2,0 x 115	14	20
6062610215	2,0 x 150	14	20
6062610220	2,0 x 200	14	20
6062600259	2,5 x 110	13	20

<b>Canulas, special with mandrin with Quincke grinding</b>			
6062570478	0,7 x 40	22	24
6062570775	0,7 x 70	22	24
6062570973	0,7 x 90	22	24
6062570980	0,7 x 120	22	24
6062570999	0,9 x 90	20	24
6062571000	0,9 x 120	20	24
6062571020	0,9 x 150	20	24
6062571120	1,2 x 100	18	24
6062572020	1,2 x 150	18	24
6062510169	0,6 x 25	23	24

# Myringotomy Kit

Several different sets for Paracentesis surgery are available

- Sterile, single use
- Customisation available

Code	Description
<b>Medi 66143300</b>	<b>Myringotomy Kit Paediatric</b>
Medi 66500140	Metal suction tip 1.4 x 70 mm
Medi 66500208	Metal suction tip 2.0 x 70 mm
Medi 66500250	Metal suction tip 2.5 x 80 mm
Medi 66500992	Suction control handle
Medi 66121370	Myringotomy scalpel
Medi 6625	Aural speculum 2.5 mm
Medi 6640	Aural speculum 4 mm
Medi 6650	Aural speculum 5 mm
Packed in a tray with 2 non-woven gauzes and a sterile field	




<b>Medi 66143400</b>	<b>Myringotomy Kit Adult</b>
Medi 66500140	Metal suction tip 1.4 x 70 mm
Medi 66500208	Metal suction tip 2.0 x 70 mm
Medi 66500400	Metal suction tip 4.0 x 80 mm
Medi 66500992	Suction control handle
Medi 66121370	Myringotomy scalpel
Medi 6640	Aural speculum 4 mm
Medi 6650	Aural speculum 5 mm
Medi 6660	Aural speculum 6 mm
Packed in a tray with 2 non-woven gauzes and a sterile field	



## Otology Suction Sets


- Ready packed Suction Sets
- Single-use, sterile, single packed in peel-pack
- Sterilisation: EO
- Material: White ABS-plastic, PVC / DEHP-free PVC, PP, medical-grade stainless steel

### Middle Ear Surgery Sets (DEHP-free PVC)



Code	Description	
6066502070	<p>Metal Suction Tip</p> <p>2 pcs 1,0 x 70 mm, brown (6066500100)</p> <p>2 pcs 1,4 x 70 mm, white (6066500140)</p> <p>2 pcs 2,0 x 80 mm, green (6066500200)</p> <p>2 m Suction tube</p> <p>Suction Tip Handle with suction control (606650992)</p>	
6066532070	<p>Metal Suction Tip</p> <p>2 pcs 1,0 x 70 mm, brown (6066500100)</p> <p>2 pcs 1,4 x 70 mm, white (6066500140)</p> <p>2 pcs 2,0 x 80 mm, green (6066500200)</p> <p>3,5 m Suction tube</p> <p>Suction Tip Handle with suction control (606650992)</p>	
6066502075	<p>Metal Suction Tip</p> <p>2 pcs 1,4 x 70 mm, white (6066500140)</p> <p>2 pcs 2,0 x 70 mm, green (6066500208)</p> <p>3,5 m Suction tube</p> <p>Suction Tip Handle with suction control (606650992)</p>	




### Mastoidectomy Suction Set (DEHP-free PVC)

Code	Description	
6066502704	Metal Suction Tip 2 pcs 1,0 x 70 mm, brown (6066500100) 2 pcs 1,4 x 70 mm, white (6066500140) 2 pcs 2,0 x 70 mm, green (6066500208) 2 pcs 3,0 x 80 mm, blue (6066500300) 3,5 m Suction tube Suction Tip Handle with suction control (606650992)	

### Stapedectomy Suction Sets (DEHP-free PVC)

Code	Description	
6066502703	Metal Suction Tip 2 pcs 0,7 x 70 mm, black (6066500070) 2 pcs 1,0 x 70 mm, brown (6066500100) 2 pcs 1,4 x 70 mm, white (6066500140) 2 pcs 2,0 x 70 mm, green (6066500208) 3,5 m Suction tube Suction Tip Handle with suction control (606650992)	
6066502280	Metal Suction Tip 2 pcs 0,7 x 70 mm, black (6066500070) 2 pcs 1,0 x 70 mm, brown (6066500100) 2 pcs 1,4 x 70 mm, white (6066500140) 2 pcs 2,0 x 80 mm, green (6066500200) 3,5 m Suction tube Suction Tip Handle without suction control (606650509)	

### Suction Sets

Code	Description	
6066502705	Metal Suction Tip 1 pc 0,7 x 70 mm, black (6066500070) 1 pc 1,0 x 70 mm, brown (6066500100) 1 pc 1,4 x 70 mm, white (6066500140) 1 pc 2,0 x 70 mm, green (6066500208) Suction Tip Handle with suction control Wullstein (606650991)	

Code	Description	
6066502706	Metal Suction Tip 1 pc 0,7 x 70 mm, black (6066500070) 1 pc 1,0 x 70 mm, brown (6066500100) 1 pc 1,4 x 70 mm, white (6066500140) 1 pc 2,0 x 70 mm, green (6066500208) 3,5 m Suction tube Suction Tip Handle with suction control Wullstein (606650991)	
6066502707	Metal Suction Tip 2 pcs 0,4 x 70 mm, lilac (6066500040) 2 pcs 0,7 x 70 mm, black (6066500070) 2 pcs 1,0 x 70 mm, brown (6066500100) 2 pcs 1,4 x 70 mm, white (6066500140) 2 pcs 2,0 x 70 mm, green (6066500208)	
6066502708	Metal Suction Tip 1 pc 0,7 x 70 mm, black (6066500070) 1 pc 1,0 x 70 mm, brown (6066500100) 1 pc 1,4 x 70 mm, white (6066500140) 1 pc 2,0 x 70 mm, green (6066500208) 1 pc 1,0 x 100 mm, brown (6066509100) 1 pc 1,4 x 100 mm, white (6066509140) Suction Tip Handle with suction control YakiSaki (606650990)	
6066502709	Metal Suction Tip 1 pc 0,4 x 70 mm, lilac (6066500040) 1 pc 0,7 x 70 mm, black (6066500070) 1 pc 1,0 x 70 mm, brown (6066500100) 1 pc 1,4 x 70 mm, white (6066500140) 1 pc 2,0 x 80 mm, green (6066500200)	
6066502710	Metal Suction Tip 1 pc 0,7 x 70 mm, black (6066500070) 1 pc 1,0 x 70 mm, brown (6066500100) 1 pc 1,2 x 70 mm, pink (6066500120) 1 pc 1,4 x 70 mm, white (6066500140) 1 pc 1,7 x 70 mm, grey (6066500170) 1 pc 2,0 x 70 mm, green (6066500208) 3,5 m Suction tube Suction Tip Handle with suction control (606650992)	

# LARYNGOLOGY

- Metal Suction Tip - Tons (Transparent)
- Metal Suction Tip - Bronch (Transparent)
- Plastic Suction Tips

## Metal Suction Tip - Bronch

- Sterile, single use
- For use in endoscopic surgery, removing foreign objects
- Transparent handle of PMMA\* for good and steady gripping
- With or without suction control



\*PolyMethylMetAcrylate

## Metal Suction Tip - Bronch

With transparent handle with suction control

Code	Size	Gauge
Medi 61616232VC	2.0 x 300 mm	14G
Medi 61616233VC	3.0 x 300 mm	11G
Medi 61616234VC	4.0 x 300 mm	8G
Medi 61616242VC	2.0 x 400 mm	14G
Medi 61616243VC	3.0 x 400 mm	11G
Medi 61616244VC	4.0 x 400 mm	8G
Medi 61616252VC	2.0 x 500 mm	14G
Medi 61616253VC	3.0 x 500 mm	11G
Medi 61616254VC	4.0 x 500 mm	8G
Medi 61616262VC	2.0 x 550 mm	14G
Medi 61616263VC	3.0 x 550 mm	11G
Medi 61616264VC	4.0 x 550 mm	8G



### With transparent handle without suction control

Medi 61616232	2.0 x 300 mm	14G
Medi 61616233	3.0 x 300 mm	11G
Medi 61616234	4.0 x 300 mm	8G

## Metal Suction Tip - Tons

With transparent handle and with suction control

Medi 66502040	4.0 x 190 mm	8G
---------------	--------------	----



## Plastic Suction Tips

Code	Size	Gauge	Colour	Description
Medi 66101023	2.0 x 100 mm	14G	Green	Single bent, short
Medi 66101023BP	2.0 x 100 mm	14G	Green	Single bent, short
Medi 66101627	2.0 x 160 mm	14G	Green	Double bent
Medi 66100627	2.0 x 160 mm	14G	Green	Single bent
Medi 66101635	3.0 x 160 mm	11G	Blue	Double bent
Medi 66100635	3.0 x 160 mm	11G	Blue	Single bent
Medi 66100635BP	3.0 x 160 mm	11G	Blue	Single bent
Medi 66101643	4.0 x 160 mm	8G	White	Double bent
Medi 66100643	4.0 x 160 mm	8G	White	Single bent
Medi 66100643BP	4.0 x 160 mm	8G	White	Single bent






# RHINOLOGY

- Metal Suction Tip - Rhino
- Metal Suction Tip - Rhino (Transparent)
- Metal Suction Tip with 2 sideholes (Transparent)
- Metal Suction Tip - Ferg (Transparent)
- Metal Suction Tip - Endo (Transparent)
- Metal Suction Tip - Fraz (Transparent)
- Metal Suction Tip - Zoellner (Transparent)
- Nasal Speculum

## Metal Suction Tip - Rhino

Made of medical grade stainless steel (virgin steel)

The nasal suction tip is shaped to a certain angle for ease of sight when used



Patient friendly  
Beveled and smooth tip end

Colour coded hub made of polypropylene that fits on all standard Luer slip parts

## Metal Suction Tip - Rhino

- Sterile, single use
- For rhinological examination, suction and flushing
- For dosing and application of salves and medicines
- 3 sizes available

Code	Size	Gauge	Colour	Description
Medi 66516020	2.0 x 150 mm	14G	Green	Luer lock
Medi 66516030	3.0 x 150 mm	11G	Blue	Round hub
Medi 66516040	4.0 x 150 mm	8G	Yellow	Round hub
Medi 66502150	2.0 x 150 mm	14G	Green	Round hub




*# The items highlighted are kept in stock, others if ordered may have a short wait time on ordering*

## Metal Suction Tip - Rhino Transparent handle

Made of medical grade stainless steel (virgin steel)

The nasal suction tip is shaped to a certain angle for ease of sight when used



Patient friendly  
Beveled and smooth tip end

Transparent handle of PMMA\*  
Thumbgrip, built in suction control  
Fits any standard universal tube end

\*PolyMethylMetAcrylate

## Metal Suction Tip - Rhino Transparent handle

- Sterile, single use
- For rhinological surgery, suction and flushing
- Transparent handle made of PMMA for good and steady gripping
- With or without suction control

### With transparent handle and suction control

Code	Size	Gauge
Medi 66509120	2.0 x 110 mm	14G
Medi 66509125	2.5 x 110 mm	12.5G
Medi 66509130	3.0 x 110 mm	11G
Medi 66509135	3.5 x 110 mm	9.5G
Medi 66516120	2.0 x 150 mm	14G
Medi 66516130	3.0 x 150 mm	11G
Medi 66516140	4.0 x 150 mm	8G



### With transparent handle without suction control

Medi 66516121	2.0 x 150 mm	14G
Medi 66516131	3.0 x 150 mm	11G
Medi 66516141	4.0 x 150 mm	8G


### With transparent handle, 2 sideholes on tip end and suction control

Medi 66509020	2.0 x 110 mm	14G
Medi 66509025	2.5 x 110 mm	12.5G
Medi 66509030	3.0 x 110 mm	11G
Medi 66509035	3.5 x 110 mm	9.5G

*# The items highlighted are kept in stock, others if ordered may have a short wait time on ordering*

## Metal Suction Tip With or without sideholes and suction control

- Sterile, single use
- For delicate procedures close to the dura, ex., or where tissue damage is to be avoided
- 2 sideholes on the tip end to avoid biopsy effect
- Transparent handle of PMMA\* for good and steady gripping
- With suction control




\*PolyMethylMetAcrylate



## Metal Suction Tip - Ferg

Made of medical grade stainless steel (virgin steel)

The nasal suction tip is shaped to a certain angle for ease of sight when used



Patient friendly  
Beveled and smooth tip end,  
available with 1 side hole

Transparent handle of PMMA  
Thumbgrip, built in suction control  
Fits any standard universal tube end

## Metal Suction Tip - Ferg

- Sterile, single use
- For rhinological surgery, suction and flushing
- Transparent handle of PMMA<sup>+</sup> for good and steady gripping
- With and without 1 side hole at the tip end
- With suction control

### With transparent handle and suction control

Code	Size	Gauge
Medi 66501020	2.0 x 150 mm	14G
Medi 66501030 *	3.0 x 150 mm	11G
Medi 66501040 *	4.0 x 150 mm	8G

### As above but with sidehole on tip end and suction control

Medi 66511020	2.0 x 150 mm	14G
Medi 66511030	3.0 x 150 mm	11G
Medi 66511040	4.0 x 150 mm	8G




\*also available in big-pack +PolyMethylMetAcrylate

# The items highlighted are kept in stock, others if ordered may have a short wait time on ordering

## Metal Suction Tip - Endo

The endonasal suction tip is bent towards the tip end for easy entrance to the sinus

Made of medical grade stainless steel (virgin steel)



Patient friendly  
Tissue saving, olive shaped tip

Transparent handle of PMMA Pengrip, no suction control  
Fits any standard universal tube end

## Metal Suction Tip - Endo

- Sterile, single use
- For sinus surgery and suction
- Specialised for FESS
- Extra softly shaped olive tip
- Non-traumatic
- Transparent handle of PMMA<sup>+</sup> for good and steady gripping
- Without suction control

### With transparent handle and without suction control

Code	Size	Gauge
Medi 66503020 *	2.0 x 150 mm	14G
Medi 66503025	2.5 x 150 mm	12.5G
Medi 66503030	3.0 x 150 mm	11G
Medi 66503035	3.5 x 150 mm	9.5G
Medi 66503040	4.0 x 150 mm	8G



\*also available in big-pack <sup>+</sup>PolyMethylMetAcrylate

## Metal Suction Tip - Fraz

- Sterile, single use
- Soft extended bend
- Transparent handle of PMMA with suction control for good and steady gripping
- Without suction control
- Longer than normal metal suction tips, i.e. more working length

### With transparent handle and suction control

Code	Size	Gauge
Medi 66508020 *	2.0 x 110 mm	14G
Medi 66508030	3.0 x 110 mm	11G
Medi 66508040	4.0 x 110 mm	8G



\* also available in big-pack

## Metal Suction Tip - Zoellner

- Sterile, single use
- The bend is close to the tip
- Transparent handle of PMMA\* for good and steady gripping
- With or without suction control

### With transparent handle

Code	Size	Gauge	Description
Medi 66504020	2.0 x 130 mm	14G	Without suction control
Medi 66504025	2.5 x 130 mm	12.5G	Without suction control
Medi 66504030	3.0 x 130 mm	11G	Without suction control
Medi 66514025	2.5 x 90 mm	12.5G	With suction control




\*PolyMethylMetAcrylate



## Nasal Speculum


- Hygienically packed, single use
- Reversed spring effect allows the physician to work with 2 hands
- White ABS plastic

Code	Description
Medi 66516100	Nasal Specula



# Preformed Sinus Catheter Post-Operative Maxillary Sinus Flusher

Code	Description
Medi 66505	Sterile Preformed Sinus Catheter



## Suction Tubes

- Sterile, single use
- Easy to connect with the suction devices and vacuum pump
- 3 versions:
  - Ready for Medioplast<sup>®</sup> Suction Handle
  - With Medioplast<sup>®</sup> Plastic Suction Tip Handle
  - With universal connectors at both ends



Length: 2.0 and 3.5 m

## ENT ACCESSORIES

### PRE CUT TUBES

Mediplast Suction Tubes, DEHP-free PVC, precut, available in 2.0 m and 3.5 m lengths, comes with female-female connectors, or female-plastic suction tip connector, or female- extension connector (for suction control handle).

Mediplast suction tubes offer a very supple and soft feel for doctors who work with long operational procedures, but still are resistant to collapse from the vacuum.

### FLUSHER SYRINGE (SJECT), 100 ML

A multi-purpose tool with an easy “one hand grip” usage. Use for flushing of wounds, for flushing ears, and also a perfect match for our Antrum Irrigation Set.

### HOSPITAL APPLICATORS

Cotton sticks, available with aluminium or stainless steel shaft. Used for applying salves and cleaning of the ear canal. Also now available in a sterile version, 5 pcs per peel-pouch, (only stain-less steel) to be used during ENT Surgery.



### TUBING ON ROLL

Mediplast can offer DEHP-free PVC tubes on roll with ID's from 4 mm to 10 mm and OD's from 6 mm to 12 mm in lengths of 30, 40 and 50 m pending on diameter size. This tube easily fed through a round top lid and cut to desired length.

### MEDIPLAST EAR FLUSHING TIP / NOZZLE

A handy tool connected to the conical end of the Flusher Syringe, used, for instance, when flushing ears, it has a plate on the cone which prohibits the cone to enter too far into the ear canal.



### VACUTIP™

A single-use, non-sterile suction device designed to safely and easily remove difficult foreign objects from the ear and nasal cavities. The suction cup conforms to almost any shape and works with a standard suction tubing and system. It fits on all Mediplast 2 mm diameter metal suction tips, enabling easy handling.



January 2016





## Endotherapeutics

Suite 1, Level 1, 51 Rawson Street, Epping NSW 2121

Phone: 1300 00 ENDO (3636) or (02) 9869 2868

Fax: (02) 9869 7723

Email: [info@endotherapeutics.com.au](mailto:info@endotherapeutics.com.au)

[www.endotherapeutics.com.au](http://www.endotherapeutics.com.au)