

o o o

ovesco
innovation in scope

LiftUp[®]

The durable solution.



Try this with other injection solutions!
www.ovesco.com



Proudly distributed by



○ ○ ○ o v e s c o
innovation in scope

LiftUp[®]. The durable solution.

LiftUp[®] creates a submucosal gel cushion at body temperature due to its thermoreversible properties. It is a new injection solution for safe and easy¹ endoscopic resection, such as polypectomy, EMR and ESD. The pronounced safety cushion lasts even after mucosa incision.²

Tissue layers and structures are well separated and exposed, thereby facilitating the resection procedure. The long-lasting lifting makes re-injection less necessary and saves time.³ Each vial contains 20 ml of LiftUp[®] pre-mixed with methylene blue and is ready-to-use.

1 Wedi, E., Ho, C.-N., Conrad, G., Weiland, T., Freidinger, S., Wehrmann, M., Meining, A., Ellenrieder, V., Gottwald, T., Schurr, MO., Hochberger, J. (2019). **Preclinical evaluation of a novel thermally sensitive co-polymer (LiftUp) for endoscopic resection.** Minim Invasive Ther Allied Technol, 28(5), 277–284. doi: 10.1080/13645706.2018.1535440

2 Wedi, E., Koehler, P., Hochberger, J., Maiss, J., Milenovic, S., Gromski, M., Ho, C.-N., Gabor, C., Baulain, U., Ellenrieder, V., Jung, C. (2019). **Endoscopic submucosal dissection with a novel high viscosity injection solution (LiftUp) in an ex vivo model: a prospective randomized study.** Endosc Int Open, 07(05). doi: 10.1055/a-0874-1844

3 Meier, B., Wannhoff, A., Klinger, C., & Caca, K. (2019). **Novel technique for endoscopic en bloc resection (EMR+) – Evaluation in a porcine model.** World J Gastroenterol, 25(28), 3764–3774. doi: 10.3748/wjg.v25.i28.3764



Further information can be found on
www.ovesco.com

Tel. 1300 00 ENDO | Fax+61 2 9869 7723
E-mail info@endotherapeutics.com.au
www.endotherapeutics.com.au

