




gastroduodenal FTRD®

Full-thickness resection now also available for the upper GI tract

- endoscopic full-thickness resection and diagnostic tissue acquisition in the stomach and duodenum
- complete set with insertion aid for the use in the upper GI tract
- transluminal and minimally invasive technique





..... gastroduodenal FTRD®


A smaller FTRD® System for endoscopic full-thickness and deep wall resection (in particular in the stomach) and diagnostic tissue acquisition in the stomach and duodenum.

For an easier and safer insertion of the system into the upper GI tract (passage of esophagus and pylorus) the gastroduodenal FTRD® set is delivered with an insertion balloon and guide wire.


The gastroduodenal FTRD® can be used for the following indications:

- not pretreated non-lifting adenoma
- non-lifting recurrent adenoma
- small subepithelial tumors
- early carcinoma


Measurements



gastroduodenal FTRD® application cap



insertion balloon



guide wire


Application aid

For a better tissue mobilization with regard to submucosal indications (especially in the stomach) a specially designed Anchor is available.



..... Application

EFTR of a GIST in the stomach with the gastroduodenal FTRD®



Source: PD Dr. A. Schmidt, Klinik für Innere Medizin II, Gastroenterologie, Hepatologie, Endokrinologie und Infektiologie, Universitätsklinikum Freiburg

Proudly distributed by



1300 00 ENDO
 info@endotx.com.au
 endotx.com.au

Ovesco Endoscopy AG
 Friedrich-Miescher-Straße 9
 72076 Tuebingen
 Germany

Phone + 49 (0)7071 96528-160
 Fax + 49 (0)7071 96528-260
 E-Mail service@ovesco.com



Ovesco Endoscopy AG is a medical device company specializing in the fields of flexible endoscopy and endoluminal surgery. Ovesco develops, manufactures, and markets innovative products for the treatment of gastrointestinal disease. Ovesco products stand for therapeutic efficacy and efficient application.

www.ovesco.com