



n&s

## Testicular Prosthesis

High technology  
and maximum safety

**Promedon**  
— People + Innovation



Testicular Prosthesis

## High Technology

The state-of-the-art silicone elastomer used to manufacture the prosthesis made the solid design resembling a natural consistency possible thanks to the low-degree hardness of the material.

## Maximum Safety

N&S can be easily implanted, causing no irritation and thus guaranteeing a safe surgical procedure, even when punctures occur. Moreover, no leakages occur due to its solid design.

### ORDERING INFORMATION

Ref.	A (mm)	B (mm)
T-Small	31	23
T-Medium	37	28
T-Large	42	32
T-XLarge	47	34

Proudly distributed by



1300 00 ENDO (3636)  
 info@endotherapeutics.com.au  
 www.endotherapeutics.com.au  
 Suite 301, 2 Banfield Road  
 Macquarie Park NSW 2113

