

UROMED UROject Injection Cannula

REF UMX 6240

Needle length 5 mm
Length of shaft 650 mm


Flexible injection cannula for toxin injection into the bladder wall

- suitable for rigid as well as for flexible cystoscopes
- reinforcement of distal material made for good guidance
- excellent toxin dosing because of empty volume and very needle
- strong haptics with easy, secure syringe attachment and removal thanks to the Luer-lock connection
- needle protection for effective insertion into the instrument
- consistently sharp bevelled needle allows accurate and effective injection

Proudly distributed by



1300 00 ENDO
info@endotx.com.au
endotx.com.au



REF	Ch./Fr.	Ø Shaft	Length of shaft	Ø Cannula	Needle length	CPN
UMX 6240	05	1,7 mm	650 mm	27 G / 0,4 mm	5 mm	12364581

Packaging: sterile
Original carton: 10 pieces